The Balanced Body

Allegro® Reformer Allegro® Tower System

THIS IS A MEDICAL DEVICE

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
Statement of Intended Use

The Allegro® Reformer Allegro® Tower System device is intended by the manufacturer to be used by people for the purposes of physiotherapy, rehabilitation, the prevention, treatment and relief of injuries, supportive, functional therapy for a disability or for general strengthening and improving the mobility of the human body. Should injury occur during usage please report to the manufacturer and competent authority of the Member State in which the user and/or patient is established.

For printed instructions contact Balanced Body.

EMERGO EUROPE
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands

IMPORTANT:

This manual is intended for medical and fitness professionals, or persons with experience in the use of this equipment. If there is a question regarding appropriateness of a particular movement, please consult a licensed health professional.

Safety Note: Warning – The Allegro contains flammable materials, please keep away from direct heat/exposed flame.

INTRODUCTION

The Allegro system is designed to combine the functions of the Allegro Reformer and the Tower or Half Trapeze into one compact and portable piece of Pilates equipment. The Allegro Tower can be installed on any Allegro and allows additional exercises to be performed without increasing the footprint of the apparatus. The following manual outlines the features of the Allegro Reformer and Allegro Tower as well as suggested exercises.

Balanced Body is the world's leading resource for Pilates equipment, education and information. The Balanced Body Studio Reformer is the biggest selling Reformer on the planet, and its Allegro Reformer now sets the Pilates standard for health and fitness clubs around the globe.

SET-UP TIPS (READ FIRST)

- » If you are setting up an Allegro Reformer with no Leg Kit, proceed with these Reformer instructions beginning on page 9.
- » If you are setting up an Allegro Reformer with with Leg Kit, install the legs first, see page 13, then set up the Reformer beginning on page 9.
- » If you are setting up an Allegro Reformer with a Tower System but no Leg Kit follow these instructions for the Reformer on page 9 and then continue to the Tower System section on page 17.
- » If you are setting up an Allegro Reformer with a Tower System and Leg Kit, install the Leg Kit first, see page 13, followed by the Tower on page 17, then the Reformer on page 9. To save time, do not re-install the head end rail cover plates as described at the end of the Leg Kit install instructions. You will need them off to install the Tower.
- » If you are setting up an Allegro Stretch Reformer, please see page 22. For Allegro Stretch footbar features please see page 25.

Features of the Balanced Body® Allegro® System

HEADREST

The headrest is used to optimally support the user's head, neck and shoulders while lying in a supine position. The Allegro has three headrest positions that are adjusted by a support block underneath the headrest:

- » **Low (flat)** – The support block is folded toward the top of the headrest. Used for clients with relatively flat thoracic spines and shallow ribcages for leg and footwork, and for any supine exercise.
- » **Safety Note:** The flat headrest position is used for all clients in exercises where they will be rolling up on to their shoulders. A flat headrest will keep the client from over flexing the cervical spine and injuring the neck.
- » **Medium** – The support block rests on the notch in the middle of the support block.
- » **High (up)** – The bottom of the support block rests on the carriage. Used for clients with a forward head or a deep rib cage to facilitate correct alignment.

Instructor Note: A towel can also be used in addition to or instead of the headrest to adjust the height of the head.

EASY TO USE PLUNGER LOCKOUT KNOBS

Balanced Body's Allegro® Reformer, and Clinical Reformer® and Studio Reformer® with Infinity footbars come with plunger lockout knobs that make adjusting and moving the Reformer footbar along the frame much easier.

To use, pull the plunger knob straight out so that the recessed cross-pin comes out of the slots as shown in Figure A. Then turn the knob so that the cross-pin rests on the bolt head as shown in Figure B. Repeat on other side.

With the knobs in the "locked out" position, reposition the footbar along the frame of your Reformer as desired. Then turn the knobs so that the cross-pin is once again recessed, firmly locking the footbar into its new position.


Figure A


Figure B

HORIZONTAL FOOTBAR ADJUSTMENTS

The Allegro is equipped with an adjustable footbar that can be moved in order to accommodate users of different heights. To move the footbar, pull the round black knobs on the bottom of the footbar out and turn 90 degrees to lock out. Slide the footbar forward or backward until it is lined up with the desired hole on the track. Turn knob until the pin engages into a hole. Make sure the pin is fully engaged.

- » **Hole 1** is closest to the footbar end of the Reformer and creates the longest distance between the shoulder rests and the footbar. It is used for taller users (over 6' or where decreased knee and hip flexion is desired).
- » **Hole 2** is for users between 5'9" and 6'
- » **Hole 3** is for users between 5'5" and 5'9"
- » **Hole 4** is for users between 4'10" and 5'5"

The heights given here are suggestions. The footbar adjustment should allow the user to have slightly less than 90° of hip flexion when the carriage is all the way in for leg and footwork.

VERTICAL FOOTBAR ADJUSTMENTS

To adjust the footbar vertically, squeeze the handle at the bottom of the footbar until the pin disengages from the plate. Move the footbar to the desired height and release the pin into the hole.

- » **High Bar** – Move the pin into the highest hole and release the handle. The high bar position shortens the space between the shoulder rests and the footbar. This position is used for foot and legwork with shorter users, and with users who have difficulty keeping their back placement due to increased lumbar lordosis or a tight back.
- » **Middle Bar** – Place the pin in the 2nd hole from the top on the plate. The middle bar position lengthens the space between the shoulder rests and the footbar. This decreases the flexion of the knees in foot and legwork and places the torso in a neutral standing position.
- » **Low Bar** – Place the pin in the 3rd hole on the plate. The low bar position increases the space between the shoulder rests and the footbar to its maximum length. This decreases the flexion of the knees in foot and legwork, decreases the flexion of the torso and hips in elephant and can be useful for taller or more flexible users.
- » **No Bar** – Place the pin in the lowest hole on the plate. This position is used to move the bar out of the way for standing exercises and for exercises where the user is lying on the box.

SPRING ADJUSTMENTS

Springs are used to adjust the resistance for different exercises. There are at least 30 different resistance settings that can be used on the Allegro. Suggested weight ranges are noted under each exercise.

CHANGING THE SPRING ATTACHMENT POINTS

On the Allegro, spring tension is adjusted by attaching different combinations of springs to the pegs under the standing platform. To safely adjust the springs make sure the carriage is in the home position.

- » **Normal position** – Place the spring on the peg directly in front of it. In this position the springs are under no tension to start with. This is designated as “B” position. There are five B positions.
- » **Pre-loaded position** – Place the spring on the peg just to the left of right, under the standing platform. This position will put the springs under a small amount of tension and will increase the resistance of the springs. This is designated as “A” position. There are six A positions.

SPRING WEIGHT

Spring resistance is indicated by the approximate number of springs suggested for a specific exercise. The spring weight indicated is a recommended starting position. Individual adjustments can be made depending on user needs and the exercise.

- » 1 spring (light): Primarily used for arm work or where the carriage is providing light support
- » 2 springs (light to moderate): Used for arm work, legwork and exercises where the carriage is providing support to the user.
- » 2 – 4 springs (moderate to heavy): Primarily used for legwork and to increase resistance for stronger users.
- » All springs: Used to maximize resistance or to stabilize the carriage for the short box abdominal series.
- » No springs: Used for added difficulty in exercises where the user needs to control the carriage (kneeling abdominals, elephant, long stretch series).

SPRING PROGRESSIONS FOR THE ALLEGRO

Please note that these spring combinations represent the usual progression for a standard new machine and may vary slightly depending on the age and specific strength of your springs. The best way to find the optimal progressions for your machine is to test it yourself.

SPRING COLOR-CODING

These color codes are standard for Balanced Body Machines

Yellow – Very light

Blue – Light

Red – Medium

Green – Heavy

ROPES AND LOOPS

Ropes are adjusted according to specific exercise demands and user size. To minimize rope adjustments, double loop straps allow the user to grip the ropes at two different lengths without having to adjust them.

The rope adjustments used in the manual are:

- » Regular loops: Resistance of loop or handle is somewhat taut on shoulder rests. With the double loop straps, the user will hold the longer loop. Standard for most exercises.
- » Short loops: Resistance of loop or handle is somewhat taut on the black pegs. With double loop handles the user holds the smaller loop. Used for rowing and some arm work exercises.
- » Very short loops: Loop or handle is approximately 5 inches shorter than the headrest. Used for kneeling arm work facing the straps such as chest expansion or thigh stretch.
- » Long loops: Loop or handle is longer than shoulder rest by a length of one cotton loop. Used for long spine stretch or for feet in the straps for users with tighter hamstrings.

RISERS

Risers can be adjusted by loosening the knob that holds the pulley in the slot on the Tower and moving the pulley up or down.

THE ALLEGRO TOWER

The Allegro Tower consists of a metal frame that attaches to the head of the Allegro Reformer. Eyebolts in the frame provide attachment points for the springs. Loops, handles or a wooden Roll-down Bar can be attached to the springs creating a wide variety of exercises. The Allegro Tower also has a Push-through Bar.

SPRING ATTACHMENT POINTS

The Allegro Tower has 24 spring attachment points creating a variety of possible exercises and a limitless amount of adjustability. Each upright has 5 eyebolts on the front and 5 on the back allowing users to perform exercises from either side of the Tower. Two eyebolts are attached to the Allegro frame to create a low position, two eyebolts are attached to the sides of the arch of the Tower and one is attached at the highest point in the center. This single eyebolt is to hold the security strap in position. Do NO attach springs to this eyebolt. Common spring attachment points are listed below and in each exercise description.

In order to adjust the tension of the springs for a particular client, move the attachment point further away from the client to make the spring heavier, move it closer to the client to make it lighter.

Low: Springs are attached from eyebolts at the bottom of the Allegro frame. These springs are used to hold the Push-through Bar in position.

Middle: Springs are attached to the 3rd eye hook from the bottom.

High: Springs are attached to the high points on either side of the arch at the top of the Tower.

ALLEGRO TOWER SPRINGS

The Allegro Tower comes with 4 sets of springs as follows:

2 Sets - Short springs

Yellow – Very light Blue – Light

2 sets - Long springs

Yellow – Very light Purple – Medium

PUSH-THROUGH BAR ADJUSTMENTS

The Push-through Bar has two possible pivot points on the Tower frame. To adjust the Push-through Bar, push the button on the top of the T-pin in order to free the pin. Pull the pins out of the uprights, move the bar to the appropriate hole and re-insert the pins. Always return the Push-through Bar to a safe height if it has been moved.

SAFETY NOTE:

It is very important that the instructor be present and spotting the client whenever the Push-through Bar is in use. It is very important that the pivot point be high enough to clear the users head when they are lying under it. If a client has an especially large head, excessive thoracic kyphosis or a large nose, the bar can hit them when they are doing exercises. The Push-through Bar should never be adjusted low enough to hit a client who is lying under it. The safety strap must always be used when the Push-through Bar is sprung from below. The safety strap must be adjusted so that the angle of the Push-through Bar, when viewed from the side, is no lower than the four or eight o'clock position and will not hit the client should their feet slip off the bar.

UPHOLSTERY CLEANING AND DISINFECTING.

You can extend the life of your upholstery by keeping it clean and free of dirt, oil and perspiration. After each use, wipe down the upholstery with a solution of mild soap and water. Then wipe it down with clean water and dry with a rag.

Safety First: A guide to proper maintenance and safe use of your Pilates equipment.

For over 35 years, Balanced Body has been introducing safety-related innovations to Pilates equipment. Many of our improvements are now industry standards, resulting in Pilates equipment that's safer today than ever before.

Safety depends on proper maintenance and safe use, in addition to the quality of the equipment. This guide was created to help you use and maintain your equipment for optimum safety. Please read it through carefully and keep for future reference. If you have any questions, give us a call. **Failure to follow these instructions may result in serious injury.**

ALL EQUIPMENT

Springs

Spring inspections are critical to maintain your equipment in safe operating condition. All Balanced Body springs should be replaced at least every two years. Certain environments and usages can shorten the expected life of the springs and you may need to replace the springs more frequently. Therefore, it is very important to inspect springs on a regular basis since worn or old springs lose resilience and may break during use. Injury may result if a spring breaks during use.

During use, do not allow springs to recoil in an uncontrolled manner. This will damage the spring and shorten its expected life.

Inspect springs for gaps and kinks (weekly or monthly, depending on frequency of use). Look for gaps and kinks between the coils when the spring is at rest. It is not unusual for the spring to have a very small gap on the tapered end (a gap is sometimes created during the manufacturing process). However, there should be no gaps in the body of the spring. If you see any gaps or kinks in the body of the spring, discontinue use and replace the springs immediately. See **Figure 1**. Additionally, corrosion anywhere on the coils will shorten the life of the spring. Discontinue using the spring immediately if you see any rust or oxidation during inspection.

Figure 1

Snaps

Inspect snaps for wear (monthly). First, verify that the snap hook is working properly. If the snap hook does not retract and return properly, discontinue using the spring immediately and replace the snap. Eyebolts can cause excessive wear on snap hooks. If the hook shows a lot of wear, discontinue using the spring immediately and call Balanced Body to replace spring or snap. See **Figure 2**.

Figure 2: Good snap: no wear on hook. Bad snap: excessive wear on hook.

EYEBOLTS, NUTS AND BOLTS

Tighten all equipment bolts and screws (monthly). Verify that all eyebolts, nuts and bolts are tight. See the section titled "How to inspect and tighten nuts and bolts."

ROPES AND STRAPS

Rope and strap wear (quarterly). Ropes should be replaced if you can see the core of the rope through the outer lining, or if the ropes are flattened. Straps should be replaced as soon as any fraying is noticed. Be sure to check the sections of rope or straps that attach to the clips and run through the pulleys.

REFORMERS

Check springbar hooks or eyebolts (quarterly). Balanced Body makes two different springbar systems:

- » Revo Springbar. Make sure springbar hooks and handle are tight.

Standard Springbar. Verify that the nuts securing the springbar hooks are tight. See section titled "How to inspect and tighten nuts and bolts."

Spring rotation (quarterly). You can prolong Reformer spring life by rotating springs of the same weight each quarter. Unhook and move to another position on the springbar. Rotating springs helps them wear more evenly.

Risers on the outside. Wood risers must be installed on the outside of the frame. Risers can loosen over time, so always make sure they are tight.

Springs hooked downward under carriage. Make sure springs are hooked in a downward position. See **Figure 3**.

Secure the carriage. When your Reformer is not in use, be sure that at least two springs secure the carriage to the springbar.

Figure 3: Springs hooked downward

Default settings. Many users have a "default setting" for Reformers. At the end of a session, the user connects a prescribed number of springs in neutral tension, sets the footbar at a pre-determined height, and sets the ropes at a specified length. This ensures that the equipment is ready for the next use, and that the carriage is secured by the springs.

Footstrap under tension in box work. When using the box and footstrap, be sure the footstrap is under tension (with snaps pulling from the top of the eyebolt) before beginning the exercise. See **Figure 4**.

Figure 4: Foot strap under tension

REFORMER WHEEL AND TRACK MAINTENANCE

Clean the tracks and wheels (weekly). For smooth carriage travel and to maintain the longevity of the wheels, we recommend that you wipe down the tracks once a week.

Disconnect the springs and clean the entire length of the tracks with a soft cloth and Balanced Body Cleaner, mild soap with water or a mild commercial cleaner such as, Fantastik® or 409®. Do not use abrasive cleansers or pads, as they can damage the anodizing on the rails. To clean the wheels, hold the cloth against the wheels while you move the carriage. If you feel a bump in the ride, dirt has adhered to the surface of the rails or wheels. Clean hair and debris out of the rails. Hair can wrap around the wheel axles and eventually build up and cause wheel failure. Use tweezers to remove hair from the wheels.

Lubrication. Never spray silicone near or inside the wheels – this can wash the lubricant out of the bearings and ruin the bearings. You can purchase dry silicone at most hardware and auto parts stores. Pulleys sometimes require lubrication to stop a squeak. Direct a very quick spray of dry silicone or Teflon spray into the pulley. "Dry" silicone does not have an oil base. Oil-based ("wet") silicone and WD40 should not be used as they attract dirt. Be careful not to over spray. You may want to remove ropes to avoid getting silicone on them.

Do not lubricate the Allegro 2 rails.

Footbar supports (quarterly). For all Balanced Body footbars with footbar support brackets, verify that the pivot screw attaching the footbar support bracket to footbar is tight, but not so tight that it prevents the support from rotating freely. For Legacy Reformers, tighten the pivot bolt to secure footbar support.

Headrest (monthly). Make sure the hinge screws and bolts on your headrest are tight.

Under the Reformer (monthly). Move Reformers and make sure you clean the floor space underneath.

Standing Platform Footbar Bumpers (wood Reformers only). If your standing platform footbar bumpers (the small plastic pieces that protect the standing platform from the footbar) are broken or damaged, please call Balanced Body to replace.

TRAPEZE TABLE (CADILLAC) & TOWERS

Cotter pins removed. These pins are located in the vertical tubes that align the canopy to the frame and should be removed as soon as installation is complete. Unremoved cotter pins can tear clothing and lacerate the skin. Use pliers to remove the pins.

Figure 5: Cotter pin before and after removal from Trap Table.

Save the pins in case you need to disassemble and reassemble the table for transportation purposes. See **Figure 5**.

Push-Through Bar (PTB) with Sliders.

The PTB moves vertically to accommodate different users and exercises. The sliders on the tubes allow for this vertical movement. Make sure these sliders are clean and easy to move. Before beginning any exercise, ensure that the sliders are properly aligned with the PTB holes and locked into position. Apply a downward force to ensure. If you notice wear on the slider knob pins, please call Balanced Body to replace. Using the PTB in this condition could cause injury.

Weekly maintenance for Push-Through Bars with Sliders. Check to ensure the pins in the plunger knobs protrude and lock into the vertical tubes. Pull both knobs out and move the sliders to a different position. Release the knobs before the next hole and continue to slide the system into position. Once over a hole the pin of the knob will automatically drop into the opening. Once each slider is in a new position apply a downward force on the PTB. The pins of the sliders should not come out of the holes. If the pins do not stay in the tube holes, the pull knobs need to be replaced. Call Balanced Body for replacement parts.

Push-Through Bar (PTB) with T-pin setting. For bottom sprung exercises, if your client's head is below the PTB, use the T-pin setting in addition to the safety strap or chain. Spotting your client is highly recommended. This is important for safety.

Push-Through Bar (PTB) control. Make sure you have enough room around the trap table to safely use the PTB without fear of hitting other people. The PTB can be dangerous if not properly

used. Only trained, experienced users should use the PTB. A spotter should always maintain control of the bar with one hand. If the user should lose control of the bar, the spotter can maintain control of it.

Correct safety strap attachment. For bottom-sprung exercises, the safety strap or chain should always secure the bar.

The safety strap or chain should wrap around the PTB and the canopy frame, not the eyebolts. The strap or chain is only as strong as the weakest link, and the frame and bar are a great deal stronger than eyebolts. **Figure 6.**

Figure 6: Safety strap holding the push-through bar at 4 o'clock. The strap is secured to the PTB and canopy frame, not the eyebolts.

Spotting your client is highly recommended. This is important for safety.

Setting the PTB for bottom-sprung exercises. For bottom-sprung exercises, the safety strap should be attached so that the angle of the push-through bar is no lower than the 4 o'clock position. This limits the range of the bar and prevents it from potentially coming into contact with the user.

Using the 4th side on the PTB along with the safety strap is highly recommended to prevent injury.

Figure 7: Safety strap holding the push-through bar at 4 o'clock. The strap is secured to the PTB and canopy frame, not the eyebolts.

CHAIRS

Dismount with control. When dismantling the chair, release the pedals slowly, with control. Don't let the pedal snap back.

Spot users. When a user is standing, sitting or lying on top of the chair, there is increased risk of falling. Standing exercises, in particular, can be unstable. Spotting users will make these exercises safer.

Hourglass spring mounts. If your chair has hourglass spring mounts and the mounts do not successfully retain the springs, replace the fiber washers (they are reddish-brown in color).

Figure 8: Fiber washer

Figure 8. If your chair is a Balanced Body Split-step Pedal Chair (Combo Chair), please consider upgrading to the Cactus Springtree).

UPHOLSTERY CLEANING & MAINTENANCE

Cleaning. You can extend the life of upholstery by keeping it clean and free of dirt, oil and perspiration. After each use, wipe down the upholstery with a solution of mild soap and water. Then wipe it down with clean water and dry with a soft towel.

Disinfecting. Equipment upholstery is coated with BeautyGard®, which offers antibacterial protection. If you want additional disinfection, Balanced Body offers Balanced Body Clean™ disinfecting solution. Use of any other solution (especially those containing essential oils) will shorten the life of some equipment and is not recommended.

HOW TO INSPECT AND TIGHTEN NUTS AND BOLTS.

Use your fingers to check nuts and bolts for tightness. If you can turn the nut or bolt with your fingers, it's too loose and should be tightened. To tighten, first tighten using your fingers. Rotate nuts and bolts clockwise to tighten. Insert a screwdriver through eyebolts to hold them steady while you tighten the nuts. Then use a small wrench to tighten the nuts further. **Figure 9.**

Figure 9: Use two fingers to tighten bolts

It is recommended to check the pins on the PTB protrude and lock into the vertical tubes appropriately. To verify their function, first pull both knobs out and start moving the sliders to a different position. Release the knobs before the next hole and continue to slide the system. Once over a hole the pin of the knob will drop into the opening. Once each slider is in a new position apply a down force on the PTB directly downwards. The pins of the sliders should not come out of the holes. If the pins do not stay in the tube holes, the pull knobs need to be replaced; call Balanced Body for replacement parts.

EQUIPMENT INSPECTION AND MAINTENANCE LOG

We suggest that you keep a maintenance log for each piece of equipment. The log should include:

1. A description of the machine including the serial number, the date and place of purchase, and the manufacturer. All of this information should appear on the invoice.
2. Date and description of all required maintenance and inspections performed.
Date and description of each repair, including name and contact information for person or company performing the repair.

MAINTENANCE SCHEDULE

All Equipment	Day	Wk.	Mo.	Qtr.
Inspect springs for gaps & kinks		✓	✓	
Inspect snaps for wear			✓	
Inspect nuts & bolts for tightness			✓	
Reformers				
Clean wheels and tracks		✓		
Inspect springbar hooks/eyebolts				✓
Rotate springs				✓
Inspect ropes/straps				✓
Inspect footbar supports				✓
Inspect One-Step springbars				✓

REPLACEMENT PARTS

To order replacement parts, or if you have any questions, please call:

U.S. and Canada: 1-800-PILATES (1-800-745-2837)
United Kingdom: 0800 014 8207
Other locations: +1 916-388-2838

Fax: 916-379-9277
Email: info@pilates.com
www.pilates.com
5909 88th Street, Sacramento, CA 95828 USA

[Click here](#) for Balanced Body Patent Data.

ASSEMBLY AND MAINTENANCE PODCASTS

View our library of assembly and maintenance videos at www.pilates.com/podcasts.

How to Assemble the Allegro® Reformer

PARTS LIST

Part number	Description	Qty
950-002	Complete footbar	1
950-020	Shoulder bolts and washers	1 kit
950-228	Shoulder rest with brackets	1 pair
616-400	Locking shoulder post	2
GEN8000	Twist lock plunger knobs	2
210-070	Soft Touch Ropes	1 pair
210-064	Double loops (Soft Touch)	1 pair
950-225	Allegro riser with pulley	2
GEN9050	5/16" Allen wrench	1
ALL0060	Allegro wrench	1
SPR9070	Red spring	3
SPR9071	Blue spring	1
SPR9241	Yellow spring	1
UPH7040	Foam Pad for Standing Platform	1
GEN9282	5/32" Allen wrench	1
GEN9472	1/4"-20 Button head screw	2
GEN9054	Split washer	2
618-008	Transport bracket preinstalled with wheel	2

**IMPORTANT PRECAUTIONS:
PLEASE READ WARNING**

To reduce the risk of serious injury, read the following important precautions before using the Allegro.

- » Read all instructions in this manual before using the Allegro. Review the Setup and Safety video included with the Allegro before using the equipment. Use the Allegro only as described in these instructions and the video.
- » It is the responsibility of the owner to ensure that all users of the Allegro are adequately informed of all precautions.
- » Use the Allegro only on a level surface. Keep hands and feet away from all moving parts. When the Allegro is not in use, leave at least two springs connected to the carriage. Keep children under the age of 12 and pets away from the Allegro at all times.
- » If you feel pain, dizziness, or shortness of breath, stop exercising immediately.
- » Before beginning any exercise program, consult your physician.

INSTALL THE SHOULDER RESTS

1. Note the two protruding studs on the bottom of the shoulder rests. Lower the shoulder rest studs into the keyholes and slide away from the carriage pad. See Figure A.

Note: The shoulder rests can be installed in two positions. One position is wider, and more comfortable for broad shoulders. Swap the left and right shoulder rests to change between standard and wide configurations.


Figure A

2. Now screw the locking shoulder posts onto the exposed bolt on each shoulder rest. See Figure B.


Figure B

Using the Locking Shoulder Posts:

Turn the posts clockwise to tighten them and lock the shoulder rests into position. Over tightening the posts will make them difficult to remove, they should be snug but not too tight. To remove the shoulder rests, turn the post counter clockwise 1-2 turns and gently slide the shoulder rests toward the carriage.

Figure C shows the storage position for the shoulder rests. Note the key hole slots at the head end of the Allegro frame.


Figure C

INSTALL THE BLACK PLUNGER KNOBS

3. Look through the large threaded hole in the trunnion and line it up with one of the holes in the frame. Screw the plunger knob into the large hole until it is completely threaded in. See Figure D.
4. Repeat on other side. Tighten both plunger knobs securely with the included open end wrench.


Figure D

To move the trunnions, pull the ball of the plunger knobs away from the frame and turn 90° to lock out. To re-engage the plunger turn the knob until it snaps into one of the holes in the side of frame.

Note: Make sure both trunnion plates are locked in the same horizontal position before moving on.

INSTALL THE FOOTBAR

5. Pick up the footbar so the seam of the footbar cover is away from the carriage pad.
6. While squeezing the footbar plunger levers, guide the footbar over the outsides of the trunnion plates. Align each footbar pin with the top hole in the trunnion plate. While continuing to support the footbar, release the footbar plunger levers to engage the trunnion plate. Be sure the pins on both sides are engaged. See Figure E.


Figure E

- Pivot the footbar slightly until the large hole at the bottom of the footbar aligns with the remaining hole in the trunnion plate.
- Slide the shoulder bolt with one washer through the large footbar hole. Then screw the shoulder bolt into the trunnion plate as tightly as possible with the large Allen (GEN9050) wrench. Repeat on the other side. See Figure E.

INSTALL THE RISERS

- Remove locking pins from the frame. Insert the risers through the holes on the top of the frame. Turn the risers until the attached pulleys point straight down the frame toward the footbar. See Figure F.
- Reinsert locking pins until the metal ring touches the frame. You may need to jiggle or slightly lift or rotate the riser to get the locking pins in all the way.


Figure F

ATTACH THE ROPES

- To attach the ropes, first engage at least one spring from the carriage to the frame to keep the carriage in the home position. Unroll the ropes and separate them.
- The Allegro now comes with Soft Touch Ropes that do not have the dog-clips to attach the loops. The clip will no longer drag on or bang into the frame! Follow the images below to install your loops. See Figure G.


Figure G

- Rest the loops over the shoulder rests. Thread the other end of each rope through a riser pulley and back into the cam cleats on the carriage to adjust the length. Be sure to go through the chrome eyestraps on both sides of the cam cleats. Always push the rope firmly down into the cleats to ensure a good grip.

INSTALL THE TRANSPORT WHEEL BRACKET

- Remove the 1/4"-20 button head screw and the split washer attached to the last hole in the frame towards the footend of the reformer, as shown in Figure H, using the 5/32" Allen key. Slide the preassembled transport bracket into the channel in the rail at the foot end of the reformer. Line up the hole in the transport bracket with the hole in the frame rail. Use the 1/4"-20 button head screw and the split washer to attach the transport bracket on to the frame using the 5/32" Allen key. See Figure H for reference. Note: The wheel goes on the inside of the bracket (towards the frame) when it is mounted.


Figure H

- Repeat this for the other side of the reformer.

STORING THE ALLEGRO REFORMER

Stacking

To see a demonstration of how to stack the Allegro, watch the Setup and Safety portion of the Allegro Introductory/Level 1 DVD. Attach the springs to keep the carriage from moving. Move the footbar to the "down" position. Remove the shoulder rests and store in the slots at the head of the frame. Put the head rest down. Remove the risers and store in the holes on the sides of the frame. Place the ropes and loops inside the Allegro frame not on the carriage. Alternate the orientation of the Allegros as you stack them. Place the head end feet onto the standing platform of the Allegro below, and so on. Stack a maximum of 5 reformers with no legs (3 with legs). Do not stack Reformer with Towers.

Standing

It is absolutely critical that the Allegro footbar be set and locked in the correct position for standing storage. Move and lock the footbar into the position closest to the foot-end of the Reformer frame. Raise the footbar to its highest position. Lift the head-end of the Allegro frame to waist height. Then raise the frame to vertical. Use your foot to steady the wheeled end of the frame as you lift to vertical. The vertical Allegro will come to rest on the footbar and the transport wheels. See Figure H.

Figure H

CLEANING & MAINTENANCE

Cleaning

Wipe the carriage pad, headrest, footbar and shoulder rests with a soft cloth and a mild, non-abrasive cleaner after each use. Keep the carriage track and wheels clean from dust and dirt. Wipe the entire track with a towel regularly. Clean the footbar trunnion track weekly. Clean the frame with a mild, non-abrasive cleaner. Keep the ropes and springs clear of dust. Cotton loops can be machine-washed. Hang to dry.

Lubrication

The foot bar trunnion C-channel can be lubricated with dry silicone spray.

MAINTENANCE

» See the included Safety First Guide.

QUESTIONS?

Call Balanced Body technical support at 1-800-745-2838 (US and Canada), or +1-916-388-2838.

How to Install the Allegro® Legs Kit

PARTS INCLUDED:

Item No.	Part Number	Description	Qty
1	618-030	Allegro Leg, A	2
2	618-031	Allegro Leg, B	2
3	618-021	Leg Nut Plate	2
4	614-020	Wheel Bracket, Allegro Ext Leg	2
5	GEN8320	Allen Key, 3/16" Ball Point	1
6	GEN7321	Screw, Button Head, 5/16-18 x 3/4"	12
7	GEN7325	Screw, 1/4-20 x 3/4"	4
8	GEN9021	Nut, Nylock, 1/4-20	4
9	GEN9282	Allen Key 5/32" Short Arm	1
10	618-135	Polyurethane Foot Pads	4
11	GEN9600	Double Sided Tape	4
12	GEN6742	Flat Head Screw 10-24 x 5/8"	4
13	GEN9856	#10 Flat Washer	4
14	GEN9059	#10 Nylon Nut	4
15	GEN9284	1/8" Allen Wrench	1
16	FAK1000	Alcohol Wipe	4

TOOLS NEEDED (NOT INCLUDED):

- » Two 1/2" or adjustable wrenches
- » One 7/16" wrench
- » 3/8" wrench or socket

PREPARE THE NEW PADS

1. Remove the legs and other materials from the packaging and lay the legs on the floor with the bottom of the leg facing up, so the "L" shape face with three holes are on the floor
2. Remove the (4) Polyurethane foot pads, and (1) Alcohol wipe from the packaging. Open one wipe and thoroughly wipe down the top side of each foot pad. This is the side opposite the ribbing. If the wipe is dirty or becomes too dry use another wipe for the remaining pads. Set the pads aside and let these dry.
3. Open an Alcohol wipe to clean the bottom of each metal foot on the reformer. If needed use a secondary wipe.
4. If you have a Tower Of Power installed, remove the ropes, loosen the 4 tower knobs, and remove the tower.
5. Using the included shorter allen wrench, remove the 4 screws holding the head end cover plates and riser pins in place. Set them aside to be re-installed later. See Figure A.



Figure A

1. Attach all springs to keep the carriage stationary. Remove ropes, shoulder rests, and risers. Place the footbar (if installed) in its lowest position and check that all pins are locked so that the footbar is secure. Then, with the help of a friend, lift the Allegro from both ends and carefully turn it upside down and set it on the floor. You may want to place a pad or other floor protection beneath the Reformer.
2. Remove the 4 rubber feet from the bottom of the Reformer frame by grasping them firmly with your hand and turning counter-clockwise. You can discard the feet unless you may want to convert back to a flat Allegro configuration.


Figure D

3. Before proceeding ensure the tops of the new pads from step 2 are completely dry and no more alcohol liquid remains on the part. If the pads are not dry, the pad and tape will not bond. Peel one side of tape off the sticker and install it on a pad. Align the 45 degree corner side of the tape with the stepped corner of the foot pad. The center hole will go over the raised center circle. See Figure B.


Figure B

4. Peel off the top side of the applied tape and install this pad assembly on the bottom of a cleaned leg. See Figure C.


Figure C

5. Use one of each: #10 screw, #10 washer, #10 nut and the provided Allen wrench, plus a socket or wrench; Align the stepped corner of the pad with the opening of the leg. The raised center will go inside the center hole in the leg. Again refer to Figure C. Tighten the screw until the end just barely sticks out the top of the nut. **DO NOT OVER TIGHTEN.**

6. Repeat steps 8-10 for all other legs.
7. Insert the leg nut plates shown in Figure D into the Reformer frame until the big hole lines up with the riser hole in the frame. When it does, drop your riser through the frame and the block to hold it in place. See Figure D and E.
8. Now maneuver the leg nut plate so that the holes in Figure E line up when you look down through the hole in the frame. This will be one of the screw holes you will use to attach the legs.
9. With the risers still in place to locate the leg nut plates, place two of the legs onto the frame rail at the head end of the Reformer, lining up the holes in the leg with the three holes in the frame. The left and right legs are not the same so keep looking for one that matches the hole pattern in the frame. Now start threading all six of the 5/16-18 screws into the frame to secure the legs in place. Start with the screws that


Figure E


Figure F

go into the block that is held in place by the risers and once they are all started, tighten them. See Figure F.

Now you can remove the risers and set them aside.
NOTE: the long end of the included 3/16" allen wrench is ball shaped to allow for more maneuverability when working in small areas. Once the screws are finger tight, use the short leg of the allen wrench to tighten the screws.

1. Install the legs onto the foot end of the frame in the same way with the 6 screws. No leg nut plates are required at the foot end.
2. Using a 7/16" wrench and the included allen wrench, install the transport wheel brackets with the 1/4-20 screws and nuts as shown in Figure G. Make sure the brackets are oriented as shown. Typically these are installed at the foot end of the Reformer as shown in Figure G.

4. Re-install the head end rail cover plates on the same side you removed them from.

NOTE: If you will be installing a Tower System next, leave the head end rail cover plates off.

5. Your legs are installed and you may now carefully lift and turn your Allegro Reformer back over onto the legs. Re-install your shoulder rests, risers, and ropes.


Figure G (Footbar not shown)

NOTE: The transport wheels can also be installed at the head end which allows you to move the Reformer with Tower more easily. Just make sure the carriage is rolled to the head end of the Reformer before lifting the foot end of the machine in this configuration.

3. Using two 1/2" or adjustable wrenches, remove the transport wheels from their current location and install them onto the brackets you just installed. Leave the old transport wheel brackets in place, as that is where a footstrap can be attached. See Figure H.


Figure H

NOTE: It is critical that the washers and wheel get installed in the same order and orientation as they were on the original brackets. Do one wheel at a time and if the wheel does not spin freely after you are done, double check against the other wheel to make sure everything is in the right order.

How to Install and Use the Standing Platform Pad

Your Allegro Reformer comes with a standing platform pad. Place the pad directly on top of the standing platform to use. This pad will make the standing platform height the same as your carriage pad. See Figure A.


Figure A

Assembling the Allegro® Tower of Power™

Prior to assembly, check the parts and components received against the following list. If you are missing any parts, contact the Balanced Body customer service department at 1-800-PILATES (1-800-745-2837) in the U.S. and Canada, or +1-916-388-2838 (international).

PREPARE YOUR ALLEGRO REFORMER

If you have purchased an Allegro 14" leg set, we strongly recommend installing the legs before installing the Tower.

Remove all ropes and risers from the Allegro. If your Allegro was purchased before February 1, 2004, remove the labels and backing on the sides. You can order a new label free of charge by calling 1-800-PILATES (1-800-745-2837 in the U.S. and Canada, or +1-916-388-2838 (international).

PARTS LIST


Description	Part Number	Qty
Tower	618-024	1
Push-through bar (PTB)	TRP0047	1
T-pin	GEN9892	2
Tower Bracket, Right	950-234	1
Tower Bracket, Left	950-235	1
Adjustable pulley and knob assembly	950-131	2
5/32" Allen key, long	212-001	1
Blue trap spring	SPR9004	2
Yellow trap spring	SPR9002	2
Long spring, purple	SPR9461	2
Long spring, yellow	SPR9006	2
Roll down bar	710-010	1
Single cotton loops	101-005	1 pair
Safety strap with carabiner	210-023	1
Tower lower spring mounts	950-236	1

IMPORTANT: READ FIRST

NOTE: If you have an original version of the Allegro with the angled corner joints (Serial Number less than 36396), follow these instructions.

- a. Using the provided 5/32" Allen wrench, loosen and remove the six screws (one on each side of the Reformer and 4 on the head end) that hold the head end of the frame on. If you have legs installed you will have to remove the leg bolts that attach to the head end of the frame. Take note of how the plastic shims between the frame members are oriented and how they fit.

- b. Install the Tower brackets as described in step number 3 below. Do not tighten the bracket screws.
- c. Re-install the head end of your Reformer frame with the plastic shims and screws removed in step a. The two shorter bolts go in the sides of the frame and the 4 longer ones go in from the head end. Get all 6 screws started, then tighten the screws until they are snug; do not over tighten them. If you have legs, reinstall those bolts also.
- d. Position the Tower brackets so they line up with the edge of the frame screws on the side of the frame. Then tighten the Tower bracket screws as described in step 4.
- e. Follow the instructions from step 6 until the end.

1. Using the provided 5/32" allen wrench, remove the plastic rail end caps and locking pins at the head end of the Allegro frame. They will need to go back on the same side, so keep track of which one goes where and keep the screws with them. See Figure A.


Figure A

2. Installing the Tower Mount Brackets:

- a. With the provided 5/32" allen wrench, loosen the 6 screws on each bracket 1-2 turns. Do not remove them, you just want some free play between the plates. See Figure B..
- b. Make sure the brackets are oriented as in Figure B (use the large grey knobs and screws as reference) and slide them into the C-channel of the frame.


Figure B

3. Re-install the plastic rail end caps on the same side you removed them from. The lanyard and locking pins are no longer needed so they can be set aside with the risers. See Figure C.


Figure C

4. Make sure both brackets are slid all the way up against the plastic rail caps. Using the provided 5/32" allen wrench, tighten the 6 screws on both Tower brackets. See Figures C and D.


Figure D

5. Install the Tower frame onto the brackets.
 - a. Orient the Tower frame so that the large threaded holes at the bottom face away from the carriage.
 - b. Lift the Tower frame over the brackets and lower it around them so it rests on the ledge at the bottom of the brackets. See Figure E.
 - c. Slide the frame toward the bracket until it is flush.
 - d. With one arm or an assistant holding the Tower frame upright and in position, start threading the 4 large knobs into the Tower frame. You may have to push on the frame to get it into position with the knob, and if one is giving you trouble, just move onto the others and come back to it.

IMPORTANT NOTE: Get all four knobs started 3-4 turns into the Tower frame before tightening any of them. It is also critical that the Tower is parallel and flush against the bracket with the knobs. You may have to push/pull/twist the leg of the Tower to make it line up.

- e. Once all four knobs are started tighten them one by one.


Figure E

6. Install the adjustable pulleys by assembling the components as shown in Figure F with the pulley on the inside of the frame. Tighten the knob when the eyebolt is horizontal as shown in Figure F. Repeat on the other side.


Figure F

7. Install the lower spring mounts in the head end of the frame with the long screw and the 5/32" allen wrench. Tighten the screw with the eyebolt oriented vertically as shown in Figure G.


Figure G

8. Rethread your ropes through the new riser pulleys.

REPOSITION THE T-PIN

For shipping purposes, the Tower is packaged with the T-pins installed on the inside of the Tower frame. See Figure H.

- The T-pin should be installed on the outside of the Tower frame for use. To remove the T-pin, push the button on the handle of the pin while holding the push through bar (PTB) securely. Pull the pin free of the PTB and Tower frame. Reinsert the pin from the outside of the Tower frame. See Figure H.


Figure H

INSTALLING MAT CONVERSION

- Disconnect all carriage springs
- Move the carriage to the head end of the frame and clear the ropes out of the way.
- Place the mat in the open space between the standing platform and the carriage.
- Remove the shoulder rests from the carriage (optional).

USING THE PRE-INSTALLED PUSH-THROUGH BAR (PTB) SAFETY STRAP

Note: Always use the safety strap for bottom-sprung exercises.

Use only with qualified instructor supervision. The PTB safety strap is used only for bottom-sprung exercises. The photo at right shows the safety strap correctly positioned to secure the push through bar. The strap includes numerous loops or "contact points" so that you can adjust the "stopping point" of the PTB if it is released during exercise.

Figure I

Note: Use the Velcro portion of the safety strap at the top of the Tower frame to secure the PTB upright when not in use.

MOVING AND STANDING AN ALLEGRO WITH TOWER

Note: Always exercise caution when moving, lifting, or storing an Allegro or an Allegro with Tower. Use proper lifting technique to avoid injury.

TO MOVE:

To move a Reformer with Tower a short distance you may be able to lift the head end and roll it on the transport wheels as if you did not have a Tower. Be aware that you will be lifting the Tower and Reformer so it will be heavier than the Reformer by itself. The Tower will also be rising quite high so you will have to be careful if you do not have high ceilings.

If moving the Reformer and Tower is not possible, remove the Tower, then move the Reformer.

- Pull the ropes back through the pulleys on the riser.
- Loosen the four gray Tower knobs and remove the Tower from the Reformer frame. If this is too heavy, the springs and push thru bar can be removed first. If you do not remove the push thru bar or roll down bar, be aware that they can swing freely.
- Now with the carriage held in the home position by at least 2 springs, lift the head end of the Reformer to roll.

Note: If you have the extended legs installed and have the transport wheels at the head end, you can leave the Tower installed and lift the foot end of the machine to roll it around. Make sure the carriage is at the head end of the frame before lifting. You will take up more room as you move around in this configuration but for some circumstances this is acceptable.

TO STAND:

- Remove the Tower frame as explained in the "To Move" section. **The Tower must be removed before standing.**
- It is absolutely critical that the Allegro footbar be set and locked in the correct position for standing storage. Move and lock the footbar into the position closest to the foot-end of the Reformer frame. Raise the footbar to its highest position. Lift the head-end of the Allegro frame to waist height. Then raise the frame to vertical. Use your foot to steady the wheeled end of the frame as you lift to vertical. The vertical Allegro will come to rest on the footbar and the transport wheels.

Note: It is absolutely critical that the Allegro footbar be set and locked in the correct position for standing storage.

How to install your Allegro® Wall Security Strap

The Allegro Wall Security Strap is designed to mount on walls with standard wood framing. If you have a different type of wall construction (brick, concrete, etc.), it may require a different fastener. Check with your local hardware store or contractor.

Parts List (included)

- » 1 adjustable belt strap
- » 1 3-1/2" eye-bolt

Tools Needed (not included)

- » Drill with 1/8" bit
- » Phillips or flathead screwdriver

INSTALL THE EYE-BOLT

1. Find a stud in the wall where you would like to stand your Allegro. Make sure you have ample space on both sides of the stud as well as above for the Allegro. Using a 1/8" drill bit, pre-drill a hole into the stud. The hole should be positioned 5' to 5 1/2' above the floor. Secure the eye-bolt into the pre-drilled hole with the screwdriver until the shaft of the eyebolt is completely within the wall.

USING THE WALL SECURITY STRAP

2. Follow your Allegro or Allegro Tower instructions for proper standing procedure.
3. Move the Allegro near the wall with the footbar pointed away from the wall (perpendicular to the wall).
4. Wrap the strap ends around one rail and connect the ends like a belt, using the adjustable clip.
5. Adjust the strap to remove any slack.

NOTE: Do not try to wrap the strap ends around the entire frame - it will not fit.

IMPORTANT: Balanced Body is not responsible for any damage or injury caused by improper wall mount installation, or the use of incorrect or improper wall mounting equipment.

If you have any questions, please call Balanced Body at 1-800-PILATES (1-800-745-2837) in the U.S. and Canada, or +1-916-388-2838 (international).

How to Assemble the Allegro® Stretch Reformer


PARTS LIST

Part number	Description	Qty
950-002	Complete footbar	1
950-020	Shoulder bolts and washers	1 kit
950-228	Shoulder rest with brackets	1 pair
616-400	Locking shoulder post	2
GEN8000	Twist lock plunger knobs	2
210-057	Soft Touch Ropes, Stretch	1 pair
210-064	Double loops (Soft Touch)	1 pair
950-225	Allegro riser with pulley	2
GEN9050	5/16" Allen wrench	1
ALL0060	Allegro wrench	1
SPR9082	Red Stretch Spring	3
SPR9080	Blue Stretch Spring	1
SPR9076	Yellow Stretch Spring	1
621-008	Oversized Standing Platform	1
UPH7040	Foam Pad - Standing Platform	1
UPH7041	Foam Pad - Oversized Standing Platform	1

IMPORTANT PRECAUTIONS: PLEASE READ WARNING

To reduce the risk of serious injury, read the following important precautions before using the Allegro.

- » Read all instructions in this manual before using the Allegro. Review the Setup and Safety video included with the Allegro before using the equipment. Use the Allegro only as described in these instructions and the video.
- » It is the responsibility of the owner to ensure that all users of the Allegro are adequately informed of all precautions.
- » Use the Allegro only on a level surface. Keep hands and feet away from all moving parts. When the Allegro is not in use, leave at least two springs connected to the carriage. Keep children under the age of 12 and pets away from the Allegro at all times.
- » If you feel pain, dizziness, or shortness of breath, stop exercising immediately.
- » Before beginning any exercise program, consult your physician.

INSTALL THE SHOULDER RESTS

The stretch reformer comes with eight different positions of shoulder rests. Shoulder rests can be adjusted based on narrow or broad shoulders. The left and right shoulder rests can be swapped to achieve these configurations.

1. Note the two protruding studs on the bottom of the shoulder rests. Lower the shoulder rest studs into the keyholes and slide away from the carriage pad. See Figure A.

Note: One position is wider, and more comfortable for broad shoulders. Swap the left and right shoulder rests to change between standard and wide configurations.


Figure A

2. Now screw the locking shoulder posts onto the exposed bolt on each shoulder rest. See Figure B.


Figure B

Using the Locking Shoulder Posts:

Turn the posts clockwise to tighten them and lock the shoulder rests into position. Over tightening the posts will make them difficult to remove, they should be snug but not too tight. To remove the shoulder rests, turn the post counter clockwise 1-2 turns and gently slide the shoulder rests toward the carriage.

Figure C shows the storage position for the shoulder rests. Note the key hole slots at the head end of the Allegro frame.


Figure C

INSTALL THE BLACK PLUNGER KNOBS

3. Look through the large threaded hole in the trunnion and line it up with one of the holes in the frame. Screw the plunger knob into the large hole until it is completely threaded in. See Figure D.


Figure D

4. Repeat on other side. Tighten both plunger knobs securely with the included open end wrench.

To move the trunnions, pull the ball of the plunger knobs away from the frame and turn 90° to lock out. To re-engage the plunger turn the knob until it snaps into one of the holes in the side of frame.

Note: Make sure both trunnion plates are locked in the same horizontal position before moving on.

INSTALL THE FOOTBAR

5. Pick up the footbar so the seam of the footbar cover is away from the carriage pad.
6. While squeezing the footbar plunger levers, guide the footbar over the outsides of the trunnion plates. Align each footbar pin with the top hole in the trunnion plate. While continuing to support the footbar, release the footbar plunger levers to engage the trunnion plate. Be sure the pins on both sides are engaged. See Figure E.
7. Pivot the footbar slightly until the large hole at the bottom of the footbar aligns with the remaining hole in the trunnion plate.


Figure E

- Slide the shoulder bolt with one washer through the large footbar hole. Then screw the shoulder bolt into the trunnion plate as tightly as possible with the large Allen (GEN9050) wrench. Repeat on the other side. See Figure E.

INSTALL THE RISERS

- Remove locking pins from the frame. Insert the risers through the holes on the top of the frame. Turn the risers until the attached pulleys point straight down the frame toward the footbar. See Figure F.


Figure F

- Reinsert locking pins until the metal ring touches the frame. You may need to jiggle or slightly lift or rotate the riser to get the locking pins in all the way.

ATTACH THE ROPES

- To attach the ropes, first engage at least one spring from the carriage to the frame to keep the carriage in the home position. Unroll the ropes and separate them.
- The Allegro now comes with Soft Touch Ropes that do not have the dog-clips to attach the loops. The clip will no longer drag on or bang into the frame! Follow the images below to install your loops. See Figure G.


Figure G

- Rest the loops over the shoulder rests. Thread the other end of each rope through a riser pulley and back into the cam cleats on the carriage to adjust the length. Be sure to go through the chrome eyestraps on both sides of the cam cleats. Always push the rope firmly down into the cleats to ensure a good grip.

14. STORING THE ALLEGRO REFORMER

Stacking

To see a demonstration of how to stack the Allegro, watch the Setup and Safety portion of the Allegro Introductory/Level 1 DVD. Attach the springs to keep the carriage from moving. Move the footbar to the "down" position. Remove the shoulder rests and store in the slots at the head of the frame. Put the head rest down. Remove the risers and store in the holes on the sides of the frame. Place the ropes and loops inside the Allegro frame not on the carriage. Alternate the orientation of the Allegros as you stack them. Place the head end feet onto the standing platform of the Allegro below, and so on. Stack a maximum of 5 reformers with no legs (3 with legs). Do not stack Reformer with Towers.

Standing

If equipped with oversize standing platform, please make sure to take it out before standing the reformer. It is absolutely critical that the Allegro footbar be set and locked in the correct position for standing storage. Move and lock the footbar into the position closest to the foot-end of the Reformer frame. Raise the footbar one down from highest position. Lift the head-end of the Allegro frame to waist height. Then raise the frame to vertical. Use your foot to steady the wheeled end of the frame as you lift to vertical. The vertical Allegro will come to rest on the footbar and the transport wheels. See Figure H.

Figure H

CLEANING & MAINTENANCE

Cleaning

Wipe the carriage pad, headrest, footbar and shoulder rests with a soft cloth and a mild, non-abrasive cleaner after each use. Keep the carriage track and wheels clean from dust and dirt. Wipe the entire track with a towel regularly. Clean the footbar trunnion track weekly. Clean the frame with a mild, non-abrasive cleaner. Keep the ropes and springs clear of dust. Cotton loops can be machine-washed. Hang to dry.

Lubrication

The foot bar trunnion C-channel can be lubricated with dry silicone spray.

MAINTENANCE

- » See the included Safety First Guide.

QUESTIONS?

Call Balanced Body technical support at 1-800-745-2838 (US and Canada), or +1-916-388-2838.

New Footbar Features for the Balanced Body® Allegro® Stretch System

HORIZONTAL FOOTBAR ADJUSTMENTS

The Allegro is equipped with an adjustable footbar that can be moved in order to accommodate users of different heights. To move the footbar, pull the round black knobs on the bottom of the footbar out and turn 90 degrees to lock out. Slide the footbar forward or backward until it is lined up with the desired hole on the track. Turn knob until the pin engages into a hole. Make sure the pin is fully engaged.

- » **Hole 1** is closest to the footbar end of the Reformer and creates the longest distance between the shoulder rests and the footbar. It is used for taller users (over 6'3" or where decreased knee and hip flexion is desired).
- » **Hole 2** is for users between 6' and 6'3"
- » **Hole 3** is for users between 5'9" and 6'
- » **Hole 4** is for users between 5'5" and 5'9"
- » **Hole 5** is for users between 4'10" and 5'5"

The heights given here are suggestions. The footbar adjustment should allow the user to have slightly less than 90° of hip flexion when the carriage is all the way in for leg and foot work.

VERTICAL FOOTBAR ADJUSTMENTS

To adjust the footbar vertically, squeeze both handles at the bottom of the footbar until the pin disengages from the plate. Move the footbar to the desired height and release the pin into the hole.

- » **Vertical Bar** – Move the pin into the highest hole and release the handle. This position is the shortest space between shoulder rests and footbar. This position is used for foot and leg work with shorter users.
- » **High Bar** – Move the pin into the 2nd hole from the top and release the handle. The high bar position shortens the space between the shoulder rests and the footbar. This position is used for foot and leg work with shorter users, and with users who have difficulty keeping their back placement due to increased lumbar lordosis or a tight back.
- » **Middle Bar** – Place the pin in the 3rd hole from the top on the plate. The middle bar position lengthens the space between the shoulder rests and the footbar. This decreases the flexion of the knees in foot and leg work and places the torso in a neutral standing position.
- » **Low Bar** – Place the pin in the 4th hole from the top on the plate. The low bar position increases the space between the shoulder rests and the footbar to its maximum length. This decreases the flexion of the knees in foot and leg work, decreases the flexion of the torso and hips in elephant and can be useful for taller or more flexible users.
- » **No Bar** – Place the pin in the lowest hole on the plate. This position is used to move the bar out of the way for standing exercises and for exercises where the user is lying on the box.

Please refer to manual for all other Allegro features.

Easy to Use Plunger Lockout Knobs

Balanced Body's Allegro® Reformer, and Clinical Reformer® and Studio Reformer® with Infinity footbars come with plunger lockout knobs that make adjusting and moving the Reformer footbar along the frame much easier.

To use, pull the plunger knob straight out so that the recessed cross-pin comes out of the slots as shown in Figure A. Then turn the knob so that the cross-pin rests on the bolt head as shown in Figure B. Repeat on other side.

With the knobs in the "locked out" position, reposition the footbar along the frame of your Reformer as desired. Then turn the knobs so that the cross-pin is once again recessed, firmly locking the footbar into its new position.


Figure A


Figure B

How to Install and Use the Oversized Standing Platform

If you ordered an Allegro Stretch Reformer, your machine comes with the "New Style" Tilting Oversized Standing Platform and Oversized Platform Pad. If you ordered a normal Allegro Reformer, refer below to determine which Platform to order.

Determine which style of Allegro Reformer you have.

- » If the rails near your trunion and standing platform have a slotted hole, follow new style instructions and order part number 950-266; Oversized Standing Platform, A1b.
- » If your reformer does not have slots in the rails, please follow the old style instructions (separate sheet) and order part number 720-063; Oversized Standing Platform, A1b.

NEW STYLE

INSTALL THE OVERSIZED STANDING PLATFORM

1. Lower the Tilting Oversized Standing Platform (O.S.P) down onto the mounted standing platform. There are two tabs on the O.S.P. which will drop into slots in the standing platform. The pins on the underside will be facing toward your carriage. See Figure B.
2. Once the O.S.P. has the tabs installed, pivot it forward and down so that the pins lower onto the rails and fall into the slots in the rails. See Figure C.
3. You can now use your O.S.P. To access and change springs, tilt the O.S.P. up and away from the carriage. When finished with the change, lower the O.S.P. back down onto the rails.


Figure B


Figure C

How to Install and Use the Standing Platform Pads for the Allegro® Stretch Reformer

Your Allegro Stretch Reformer comes with two standing platform pads, one for the fixed standing platform, the other for the oversized standing platform. Place them directly on top of the standing platform to use. These pads will make the standing platform height the same as your carriage pad. See Figures A and B.

Note: You must remove the smaller pad to install the O.S.P. ONLY use one pad at a time.


Figure A - Standing Platform Pad


Figure B - Oversized Standing Platform Pad

SAFETY & MAINTENANCE

SAFETY

To reduce the risk of injury or damage, read all instructions and the following important precautions before using the Allegro Reformer Allegro Tower System.

- » It is the responsibility of the owner to ensure that all users of the Allegro Reformer Allegro Tower System are adequately informed of all precautions.
- » Use the Allegro Reformer Allegro Tower System only on a level surface.
- » Keep hands and feet away from all moving parts. When the Allegro Reformer Allegro Tower System is not in use, leave at least two springs connected to the frame.
- » Keep children under the age of 12 and pets away from the Allegro Reformer Allegro Tower System at all times.
- » If you feel pain, dizziness, or shortness of breath, stop exercising immediately.
- » Do not stand on the footbar.
- » Before beginning any exercise program, consult your physician.
- » Do not stand the Allegro Reformer Allegro Tower System on end without the wheel kit.

MAINTENANCE

- » Depending on frequency of use, the springs should be safe to use for one to two years (3000 hours). Replace any spring that is kinked, bent or shows separations while at rest.
- » If your ropes need to be replaced refer to the instructions on the underside of the carriage. The Reformer will need to be turned on its side (please place a Pilates mat on the floor next to the Reformer to protect the finish), or stood on end if you have transport wheels.
- » Remove hair and debris from the axles of the wheels. Debris can create pressure on the wheels, causing them to wear. If you hear noise from the bearings, replace them.
- » Repair any tears in the upholstery with an upholstery repair kit available at auto parts stores, or have the repair done locally.
- » Contact Balanced Body for other upholstery and re-upholstery options. Replacement upholstery is available from Balanced Body.
- » Check and tighten all screws and bolts. Replace any missing screws, retention pins, or other parts.
- » Check the footbar bolts on a monthly basis to make sure they are still secure.

CLEANING

Wipe the carriage pad, headrest, footbar and shoulder rests with a soft cloth and a mild, non-abrasive mixture of soap and water after each use. Keep the carriage track and wheels free from dust and dirt. Clean the frame with a mild, non-abrasive mixture of soap and water. Keep the ropes and springs clear of dust. Cotton loops can be machine-washed. Hang to dry.

NOTE: DO NOT LUBRICATE OR GREASE ANY COMPONENTS ON THE FOOTBAR, CARRIAGE, WHEELS OR RAILS. If you have any issues or questions please contact Balanced Body Technical Support. By adding lubrication or grease it may not lock the components securely and may become unsafe.

Questions? Call Balanced Body Technical Support at 1-800-PILATES or +1-916-388-2838 or info@pilates.com.

Reformer Exercises

FOOTWORK, ALL LEVELS

10 reps

Springs: 2 – 4 springs

Bar: Middle or High

Head rest: Up

FOCUS

- » **Breathing** – exhale out/inhale in, or inhale out/exhale in
- » Neutral spine
- » Pelvic stability
- » Hip, leg and ankle alignment
- » Hip, leg and ankle strengthening
- » **Circulation**
- » **Isolation** – release unnecessary tension in the upper body and hips

PRECAUTIONS

Sensitive to ankle, knee, hip flexion, spinal compression

STARTING POSITION

Supine on carriage, feet on foot bar, legs hip width apart

Heels

Heels on foot bar, push back and return

Toes

Ball of foot on bar, heels slightly raised, push back and return

Prehensile

Ball of foot wrapped around the bar, push back and return

Pilates V

Ball of foot on bar, turned out, low releve, heels together, push back and return

Flex/Releve

Ball of foot on bar, parallel, legs straight, plantar flex ankle, dorsiflex ankle, plantar flex ankle, bend knees, push back to starting position

2nd position

Heels at ends of bar, slight turn out, push back and return

Heels


Toes

Prehensile

Pilates V

Running in place

Ball of foot on bar, dorsiflex one heel, bend the other knee, alternate legs 20-50 times

ROLLBACK, LEVEL 1

10 reps

Springs: 1 – 2

Bar: None

Loops: Short

FOCUS

- » Breath - exhale roll down/inhale roll up or inhale down/exhale up
- » Spinal flexibility
- » Abdominal strength
- » Shoulders down
- » Neck long
- » Soft hip flexors

PRECAUTIONS

Back injuries, neck injuries, hip flexor injuries. Be cautious with osteoporosis

PREREQUISITES

Comfortable in spinal flexion

STARTING POSITION

Seated facing ropes, knees bent, loops in hands with arms straight and elbows soft

Standard Exercise

Hold loops with arms straight, roll down, curving back and staying lifted, roll back up maintaining slight flexion in spine

Oblique Variation

Rotate knees to one side and torso to the other, roll down, switch rotation and roll up

Arm work Variations

Roll back and with low back on carriage add bicep curls or deltoid lifts for 3 reps, roll up

HUNDRED, ALL LEVELS

Sets: 10

Springs: 1 – 3

Bar: None

Loops: Regular

Headrest: Up

FOCUS

- » Percussive breathing - in for 5, out for 5
- » Stable pelvis – imprinted or neutral
- » Hollow abdominals
- » Abdominal strength
- » Shoulders down
- » Neck long

PRECAUTIONS

Back injuries, neck injuries, hip flexor injuries, osteoporosis

PREREQUISITES

Hundred on the mat

STARTING POSITION

Lying supine on carriage, knees at 90 degrees, hands in loops, arms to ceiling

Level 1


Knees bent at 90 degrees, reach arms to sides as the head and upper body lift off the carriage, pulse arms with breath

Level 2

Legs straight up to ceiling, reach arms to sides as the head and upper body lift off the carriage, pulse arms with breath

Level 3

Reach arms to sides as the head and upper body lift off the carriage, straighten legs to ceiling then lower legs keeping low back on mat, pulse arms with breath

COORDINATION, LEVEL 2

6 reps

Springs: 1 – 2

Bar: None

Loops: Regular

Headrest: Up

FOCUS

- » Breath – inhale start, exhale reach, inhale open/close, exhale return
- » Stable pelvis – imprinted or neutral
- » Hollow abdominals
- » Abdominal strength
- » Adductor strengthening
- » Shoulders down
- » Neck long

PRECAUTIONS

Back injuries, neck injuries, hip flexor injuries, be cautious with osteoporosis.

PREREQUISITES

Mat or Reformer Hundred

STARTING POSITION

Lie supine on carriage, knees at 90 degrees, hands in loops with elbows bent at 90 degrees and upper arms on carriage

Standard exercise

Reach arms to sides as the head and upper body lift off the carriage to Hundred position, open and close legs, bend knees first, then lower head and return arms to starting position

REVERSE ABDOMINALS, LEVEL 4

10 reps

Springs: 1 – 2

Bar: None

Loops: Regular

FOCUS

- » Breath - exhale to pull the knees in, inhale to release
- » Abdominal strength
- » Iliopsoas strength
- » Iliopsoas and abdominal coordination
- » Imprinted spine
- » Shoulders down
- » Neck long

PRECAUTIONS

Back injuries, neck injuries, hip flexor injuries, osteoporosis

PREREQUISITES

Mat Roll-up, strong abdominals, ability to maintain imprinted spine

STARTING POSITION

Supine facing ropes with an imprinted spine, knees bent at 90 degrees, loops around the knees, head supported with the hands, elbows wide

Standard Exercise

Maintaining imprinted spine, hollow out abdominals and pull knees above 90 degrees while flexing the torso and lifting the head up, return legs to starting position maintaining an imprinted spine

Oblique Variation

Rotate torso to the right, reach left hand across toward right knee, pull knees in toward the chest as arm reaches across (Repeat 4 – 8 times on one side and switch)

KNEELING ABDOMINALS, LEVELS 1-2

10 reps

Springs: 0 – 2

Bar: None

Loops: None

FOCUS

- » Breath – Exhale as knees pull in/ inhale to return to starting position
- » Pelvic stability
- » Abdominal strength
- » Pelvic stability
- » Scapular stability

PRECAUTION

Shoulder, elbow and wrist problems, inability to kneel

PREREQUISITES

Mat all fours hollowing (pregnant cat abdominals)

STARTING POSITION

Kneel on all fours facing straps with hands on the edges of the frame

Standard Exercise (flat back)

Face straps with knees against shoulder rests, keeping spine neutral and parallel to the ground, pull carriage toward the head by flexing at the hips

Variation – Oblique 1

Place both hands on one frame allowing torso to side bend and rotate while keeping hips as square as possible, pull both legs in increasing the movement of the torso

Variation – Single leg 1

Lift one leg off carriage and extend it in a straight line behind you, keeping the shoulders as square as possible pull one leg in maintaining a neutral back position

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Loops: Regular

FOCUS

- » Breath - exhale out/inhale in or inhale out/exhale in
- » Spine to mat or neutral spine
- » Hollow abdominals
- » Abdominal strength
- » Pelvic Stability
- » Hamstring, adductor and gluteal strength
- » Hamstring and adductor flexibility
- » Leg and hip alignment
- » Hip range of motion

PRECAUTIONS

Hip flexor injury, limited hamstring flexibility, back injuries, weak abdominals

PREREQUISITES

Adequate hamstring flexibility, ability to stabilize the back, Hundred

STARTING POSITION

Supine on carriage, loops around arches

Leg lowers

Hips stable, begin with inner thighs together, lower and raise legs

Variations: Legs parallel, turned out or turned in Magic circle or ball between the legs


Scissors

Hips stable, begin with inner thighs together, open legs to sides and return

Variations: Legs stay over hips, legs move down toward the bar while carriage moves (V's)

Circles

Hips stable, begin with inner thighs together, moving legs down and around in circles or D's, reverse directions Variations: Legs parallel, turned out or turned in, or knees in straps

PULLING STRAPS, LEVEL 1

4 – 8 reps

Springs: 1 - 2 springs

Box: Long

Loops: Regular

Bar: Low or none

FOCUS

- » Breath – inhale pull, exhale release
- » Hip bones to mat
- » Hollow abdominals
- » Triceps, latissimus, lower trap and upper back strength
- » Head in line
- » Scapula placement
- » Leg and torso alignment

PRECAUTIONS

Shoulder, elbow and wrist injuries, some back problems

PREREQUISITES

Ability to lie prone

STARTING POSITION

Lie prone on long box with chest off front edge of box, foot bar down, facing ropes, grasp ropes.

Arms down

Pull ropes to hips along edge of carriage

Arms out to sides

Pull ropes to hips with arms extended out to the side

Upper back lift

Pull ropes along edge of carriage and lift upper back

Triceps press


Pull ropes to hips along edge of carriage, keeping elbows in place, flex and extend the elbow

SWAN (LEVEL 2)


STARTING POSITION

Prone on long box, facing low or middle footbar, hands on bar.

Press out to straight arms, press down on the bar to lift into back extension, press torso back to flat, bend elbows.

Pulling Straps

Pulling Straps with Upper Back Lift

Swan

SHORT BOX - ABDOMINALS, ALL LEVELS

6 reps

Springs: 4 - 5

Box: Short

Strap: Foot strap

Pole in hands

FOCUS

- » Breath – inhale down/exhale back or exhale down/inhale back
- » Abdominal strengthening
- » Abdominal hollowing
- » Back strengthening
- » Back flexibility

PRECAUTIONS

Back, neck and shoulder injuries, osteoporosis

PREREQUISITES

Reformer Rollbacks, Cadillac Rollbacks

STARTING POSITION

Sitting on short box, facing footbar, feet under foot strap, knees slightly bent

Standard Exercise

Roll down with a long curve and roll back up

Oblique variation

Roll down with a long curve, rotate the torso to each side, roll back up

PELVIC LIFT, LEVEL 2

6 reps

Springs: 2 - 3

Headrest: Down

FOCUS

- » Breath – exhale up, inhale back, exhale in, inhale roll down
- » Hamstring and gluteal strengthening
- » Abdominal hollowing
- » Pelvic stability
- » Spinal mobility

PRECAUTIONS

Back problems, limited knee flexion, limit roll up with cervical problems

PREREQUISITES

Mat Pelvic Press

STARTING POSITION

Lying supine, ball of feet or heels on bar, legs parallel or V-feet and hips externally rotated

Standard Exercise

Roll pelvis up off the carriage, push back, keeping hips at one level, bring carriage in, roll down

Leg positions

Legs parallel, heels under sit bones

Legs parallel and inner thighs together

Legs in Pilates V with heels together

Heel frame width apart with the hips turned out.

ARM WORK, LEVEL 1 - 2

4 - 10 reps

Springs: 1-2

Box: long, short or none

Loops: very short, short or regular

Focus

- » Breath – inhale pull/exhale release
- » Biceps, triceps, pectoralis and deltoid strengthening
- » Scapular stabilization
- » Torso stabilization
- » Sitting posture

Precautions

Wrist, arm or shoulder problems, back problems with limited sitting ability

Sitting variations for all exercises

Sitting on carriage, cross-legged, legs straight or kneeling

Exercises Facing the Straps

Biceps

Holding very short loops in hands, bend elbows to pull straps to shoulders.

Triceps/Posterior Deltoid

Loops in hands, arms straight, pull straps back level with hips and pulse arms back

Exercises Facing Footbar

Serve a tray

Sitting facing footbar, regular loops in hands, elbows bent, reach forward, straighten arms, open arms to the side palm up and return

Hug a tree

Sitting facing footbar, regular loops in hands, arms out to sides, soft elbows, bring fingertips toward each other

Biceps

Serve a tray

Hug a tree

Triceps

STANDING, LEVEL 1 - 3

8 reps

Springs: 1 spring

Standing platform

FOCUS

- » Breath – exhale out, inhale in
- » Adductor and abductor strengthening
- » Abdominal hollowing
- » Standing alignment
- » Balance

PRECAUTIONS

Balance problems

Prerequisites

None

STARTING POSITION

Standing with one foot on standing platform, and one on the carriage. Face side of the room parallel to the side of the Allegro. Put the first foot on the frame and the second foot on the carriage.

Standard Exercise

Stretch legs open, close legs

Leg Position Variations

Legs straight, parallel, turned out

Legs bent, parallel, turned out

Allegro® Tower Exercises

DEFINITION OF EXERCISE SET-UP TERMS

- » **Level:** The level of expertise needed to undertake exercise.
- » **Reps:** How many times the exercise is performed.
- » **Springs:** Number and location of springs on the Allegro Tower.
- » **Loops:** Which loops should be used during exercise.*
- » **Focus:** What should be emphasized during exercise.
- » **Precautions:** Physical conditions that may limit or exclude a participant. Exercises may need to be modified for people with these conditions.
- » **Prerequisites:** Specific exercises that must be mastered before undertaking a new exercise.
- » **Starting Position:** Where to begin the exercise on the Allegro Tower.
*If applicable

ROLLBACKS, LEVEL 1

Springs: 2 long yellow or 2 short yellow springs from high position

Reps: 6-10

Rollback Bar or handles

Focus

- » Breathing – exhale roll down, inhale at the bottom, exhale to roll up
- » Balance between abdominals and lumbar extensors
- » Soft neck and shoulders
- » Maintain C-curve
- » Soft hip flexors

Precautions

Shoulder and neck problems, some low back problems, osteoporosis

Starting position

Sit facing Tower, holding on to bar or handles, knees soft, feet on metal bars.

Standard Exercise

Hold bar with arms straight, roll down, curving back and staying lifted, roll back up maintaining slight flexion in spine

Oblique Variations

- » Wooden bar (Water skiing)
- » Sit diagonally on table, place left foot against pole, cross right foot over ankle, place left hand on bar and reach right arm open while rotating torso to the right.
- » Roll down and up maintaining rotation, then switch sides.

Roll Backs

Water skiing

FEET IN STRAPS, LEVEL 1

Reps: 6

Set up: 2 springs

Bar: Any

Loops: Regular

Focus

- » Breath - exhale out/inhale in or inhale out/exhale in
- » Spine to mat or neutral spine
- » Hollow abdominals
- » Abdominal strength
- » Pelvic Stability
- » Hamstring, adductor and gluteal strength
- » Hamstring and adductor flexibility
- » Leg and hip alignment
- » Hip range of motion

Precautions

Hip flexor injury, limited hamstring flexibility, back injuries, weak abdominals

Prerequisites

Adequate hamstring flexibility, ability to stabilize the back, Hundred

Starting Position

Supine on carriage, loops around arches

Leg lowers

Hips stable, begin with inner thighs together, lower and raise legs

Variations: Legs parallel, turned out or turned in Magic circle or ball between the legs

Scissors

Hips stable, begin with inner thighs together, open legs to sides and return

Variations: Legs stay over hips, legs move down toward the bar while carriage moves (V's)

Circles

Hips stable, begin with inner thighs together, moving legs down and around in circles or D's, reverse directions Variations: Legs parallel, turned out or turned in, or knees in straps.

Starting Positions - Leg Lowers 1

Leg Lowers 2

Scissors 1 / Circles 1

Scissors 2

Circles 2

Circles 3

FEET IN STRAPS SIDELYING – LEVEL 1-3

6-10 reps

Springs: Long springs from middle position

Loops around arches

FOCUS

- » Breath - exhale down/inhale up
- » Correct side-lying position
(waist up, hips and shoulders in line)
- » Hollow abdominals
- » Pelvic stability and isolation of the leg from the pelvis
- » Adductor, medial hamstring and external rotation strengthening

PRECAUTIONS

Some back injuries, knee injuries, and unstable sacroiliac joints

STARTING POSITION

Lie on your side on Allegro Tower with back of body in line with back edge of mat and legs slightly forward.

Support body by bracing the bottom arm against upright pole or resting head on arm.

Place loop around arch

Standard Exercise

Pull top leg down toward bottom leg.

Maintain correct side-lying position.

Variations: parallel, turned-out (larger range of motion), turned-in.

Ovals

Maintaining correct alignment, move top leg in a small circle in both directions.

Variations: parallel, turned-out, turned-in


Front-Back Kick

Maintaining correct alignment, swing top leg forward and back (as in the Side Kick on the mat).

Variations: parallel, turned-out, turned-in

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FOOT AND LEGWORK, LEVEL 1

10 reps

Springs: 2 long purple springs from the bottom on Push-through Bar

Safety strap on

FOCUS

Breath – inhale push, exhale return

Spine to mat or neutral spine

Leg, ankle and foot alignment

Calf and hamstring flexibility

Foot, ankle and lower leg strength

PRECAUTIONS

Back injuries, knee injuries

STARTING POSITION

Supine on Allegro Tower with feet on Push-through Bar, and springs attached from low position onto the Push-through Bar

Plies

Lie supine with the Push-through Bar in line with anterior hip crease, flex knees, flex hips, with metatarsals or heels on the bar, and straighten legs and return.

Foot position variations: Heels, Toes

Leg variations: parallel, turned out, v-feet, wide 2nd position, single leg

Plie/Releve

Lie supine with the Push-through Bar in line with anterior hip crease, hips flexed, knees bent, metatarsals or toes on the bar.

Push the bar up toward the ceiling straightening the knee, plantarflex the ankle, dorsiflex the ankle and return.

Variations: parallel, turned out, single leg

Plantarflexion

Lie supine with Push-through Bar in line with anterior hip crease, legs straight, metatarsals or toes are on bar, plantarflex and dorsiflex the ankles.

Variations: parallel, turned out, single leg, running in place

Plies

Plantarflexion - 2 legs

Plantarflexion - alternating legs

SEATED PUSH THROUGH, LEVEL 1

6 reps

Springs: 1 short red spring from above on Push-through Bar

FOCUS

- » Breath – exhale stretch forward, inhale stretch up or reverse
- » Spinal flexion and extension
- » Hamstring flexibility
- » Scapula stability/mobility
- » Abdominal hollowing and lift
- » Coordination of breath with spinal mobility

PRECAUTIONS

Some back injuries, shoulder injuries

PREREQUISITES

Mat – Spine Stretch

STARTING POSITION

Sit on table facing Tower, with feet against upright bars, and both hands on the Push-through Bar (knees can be bent or straight depending on flexibility).

Standard Exercise

Push bar down, curve spine forward beginning with top of head, reach the bar forward and stretch.

Return by hollowing out abdominals and stacking one vertebra on top of another until sitting up on the sit bones.

Press bar up and lean forward from hips with a flat back before beginning again

CAT, LEVEL 3

4 reps

Springs: 2 short springs from above on Push-through Bar

FOCUS

- » Breath – exhale roll down, inhale to extend out, exhale pull back, inhale uncurl.
- » Spinal flexibility
- » Scapula stability/mobility
- » Abdominal hollowing and lift
- » Coordination of breath with full spinal mobility

PRECAUTIONS

Back injuries, shoulder injuries, knee injuries, and osteoporosis

PREREQUISITES

Mat – Cat/camel stretch

STARTING POSITION

Kneeling on table, hands on Push-through Bar with the bar close to the body.

Standard Exercise

Press bar down and roll the spine down beginning with top of head.

Reach bar away as spine elongates into extension.

Return by pulling abdominals in and curling spine back into flexion before stacking vertebra one on top of the other to return to the starting position. (Keep the hips pressed forward over the knees as much as possible.)

CIRCLE SAW, LEVEL 2

4 reps

Springs: 2 short yellow springs from above on Push-through Bar

FOCUS

- » Breath – exhale stretch, inhale to circle out, exhale reach, inhale circle back
- » Spinal flexibility in rotation
- » Scapula stability/mobility
- » Abdominal hollowing and lift
- » Coordination of breath with full spinal mobility

PRECAUTIONS

Some back injuries, shoulder injuries

PREREQUISITES

Mat - Saw

STARTING POSITION

Sit on table facing Tower, with feet against upright bars, right hand pushing up on the Push-through Bar, left hand reaching toward right foot.

Standard Exercise

Sweep left hand toward left foot and continue to circle arm out and around, allowing torso to lean back while maintaining abdominal lift and torso integrity until the left hand reaches over the right arm toward the right foot.

Reverse the circle

Do four repetitions using each arm, keep both sit bones anchored for pelvic stability (or release one slightly to increase stretch).

SEATED DOUBLE ARMS, LEVEL 1-2

4-10 reps

Springs: 1 - 2 short red springs from high position on Push-through Bar
Push-through Bar
Dowel or padded stick

FOCUS

- » Breath – exhale pull, inhale release or inhale pull, exhale release
- » Scapular stability and mobility
- » Shoulders down
- » Abdominal hollowing
- » Balance on sit bones

PRECAUTIONS

Shoulder injuries, modify for tight adductors

STARTING POSITION

Straddle table facing Tower, hands on Push-through bar, or sit cross-legged inside frame of Push-through Bar. Can also stand on outside of Allegro facing Tower.

Pull downs with stick

Hold the dowel or padded stick as wide as possible with palms facing you.

Pull bar down, focusing on soft neck and upper shoulders.

Return bar to starting position without raising the shoulders.

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