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For stiff shoulders and neck

Swedish Posture®

EN Elegant shoulder wrap filled with flaxseed and lavender for a relaxing feeling. Embrace Shoulder Wrap is comfortably used over the shoulders, upper back and neck for a warm and calming sensation, like getting a comforting hug. Heat the wrap in the microwave or chill it in the freezer for thermo therapy. The wrap can be used when feeling cold or for neck pain, shoulder stiffness, tired muscles, tension, impingement, swelling, stress and headache.

Hot therapy

- Place the weight wrap in a clean microwave and heat it for 30-60 seconds (600 W)
- Shake the flaxseed around inside the wrap.
- Repeat until the wrap is warm enough.
- The wrap should feel warm, not too hot to touch.
- The wrap MUST be supervised when heated in the microwave.
- The wrap keeps warm for about 30 minutes. Then it can be reheated.

Cold therapy

- Put the wrap in the freezer for at least 30 minutes.
- The cooling effect can reduce swelling

Material

Textile: 70% cotton, 30% Linen

Filling: Flaxseed, Lavender

Made in China

- Use when feeling cold or for neck pain, shoulder stiffness, tired muscles, tension, impingement, swelling, stress and headache
- Use hot or cold for thermo therapy
- Lavender scent gives a calming effect
- The heavy weight gives a comfortable and relaxing sensation, like a loving hug
- Use anytime when feeling tired, tensed, stressed or in pain

What is thermo therapy?

Thermo therapy involves applying heat or cold to joints/muscles. Heat improves blood circulation and relaxes muscles, while cold can numb pain, decrease swelling and constrict blood vessels to relief pain.

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Inspirera & Co AB/