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Manual instruction

Product-Description

TheraBeans are an easy to use and highly efficient exercise-device for all the muscles of feet, lower leg, hands and forearms.

Developed from a training-method of high performance sports - the "corn-box" of Russian triple-jumpers, a sport with an extreme stress for the upper joint of the ankle and the stabilizing muscles of the foot - the TheraBeans now offer a hygienically perfect utilization within the therapy of paretic muscles in arms and legs.

Every muscle is trained with a simple set of exercises without being overloaded. At the same time a massage-effect with an increase of blood-circulation is achieved. Also strong sensorimotoric responses are created - the TheraBeans therefore are an valuable tool for psychomotorical therapy too.

Priority of use is the rehabilitation of paretic muscles in upper and lower extremity (orthopedics, surgery, sport-rehabilitation, neurology) for a simple, efficient training of patients under own regiment.

Features

TheraBeans are made of highly stable, PVC-free plastic. The used colors are cadmium-free and food-proof.

The two sizes are ideally tuned for the training and offer an even, well-dosed resistance to movements.

Resistance is generated by the mechanical friction of the beans when moved and is nearly neutral in direction - thus having no load-free phase - training in the beans is therefore always very easy in the beginning, though very exhausting for the muscles in a short time due to missing relaxation periods.

For an effective training the following amounts are needed:

Base-Set foot/forearm training:
• Quantity of TheraBeans approx. 50 liters/30kg
• Case-size approx. 45cm x 40 cm x 40 cm

Training both feet/both forearms:

- Quantity of TheraBeans approx. 80 liters/50kg
- Case-size approx. 60 cm x 45 cm x 40 cm

The used case must be so large, that the TheraBeans and the body-part to be exercised can be taken into the volume even when moving.

TheraBeans are easily cleaned and disinfected. Please take care for sufficient rinsing with clear water after disinfection, to avoid skin irritations when working later.

Regular Use

TheraBeans are a training-device for a therapy of paretic muscles at various indications.

They are used in a suitable case, which is firm enough to cope with the weight of the beans and the load of the exercising person.

It must be large enough to provide ample space for the exercising person and permitting the exercise and the movement of the beans.

The case must be placed secured and slip-free. The exercising person must be able to maintain a comfortable and safe exercise-position.

Please take care for proper back-safe working when moving the filled case or use suitable transport-tools.

Preparations for the Exercises

Take off all jewelry and accessories (as possible), if needed cover irritated or wounded parts of skin.

Put arm or leg/foot, according to the desired exercise, slowly into the beans and work yourself downwards into the beans with small movements until your arm is covered with beans shortly above the elbow or your leg at least one hand below the knee.

Follow the adjoining exercise-instruction poster or the correspondent instructions of your therapist for the specific exercises. Take care for minute execution of the instructions as there is a risk to be overloaded or wrongly loaded (problem e.g. at fractures).

Always start the exercises slowly and get accustomed to the movements and the resistance (warm up).

Never work into the pain, always remain within your personal range of tolerance - dangerous overload is thus avoided reliably.

Slowly increase the intensity of the exercise (velocity/range of motion) to your personal maximum possible strength and work out until you are no longer able to execute the movement properly or your muscles are exhausted.

Finish the training and take arm/leg out of the beans. Do not take up the work-out before finishing a relevant break (approx. 5 - 8 h).

Warnings

- Children below the age of 3 are not allowed to work with and in the TheraBeans, if not monitored
- TheraBeans are not suitable to be eaten
- TheraBeans should only be used when properly instructed and on suggestion of a therapist
- Do not use TheraBeans at strong circulatory disorders (arterial occlusion etc.) and known thrombosis (risk of embolism).
- Wounded or irritated skin parts must be covered and protected
- Exact execution of the exercise movements is essential at existing fractures or other instabilities in joints and bones to avoid dangerous overload. Please consult your therapist on this topic.

Technical Data

Material: Polypropylene
Density: 0.6 kg per liter

Cleaning

Suitable for standard detergents/disinfectants for regular cleaning (sometimes surface of beans will get dim after a time).

Conditions of storage, transport and use

Ideal conditions are temperatures above 0° Celsius. Avoid direct sun. The TheraBeans contain no uv-stabilizer as this could only be achieved by recruiting plastic-softeners.

Temperature range is from 1°Celsius up to 50°Celsius, warmed up in oven or on heating (do not use microwave as the colour contains small metallic quantities).

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Jakobs GmbH