Fußwalzen

Training der Stand- und Gehsicherheit
Stability and Walking safety will be improved

Art.-Nr.: 653 098

Deutsch

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English

The Footboards moveable fulcrum can be used individually and allow varied exercises for foot-leg-axle training (foot- > ankle- > knee- > hip-joint) to improve coordination and stabilization, strengthening and stretching of the foot to erect the horizontally foot arch and train the longitudinally arch for a healthy foot. Your stability and walking safety will be improved noticeably. Proved in many cases in prevention of fall prevention. Ideal for strenghtening exercises in front leaning rest on ground when hand or forearm rests on the footboards. Also ideally used as supplemental parts to Balance Boards (fig. 2), Spring Boards or Rotating Balance Board to increase level of difficulty. Stowed away quickly in every bag, the Footboards moveable fulcrum are an optimal appliance to each training.

Item no. 653 098

Specification:

2 pcs Footboards moveable fulcrum in solid ash, 35x10x5 cm, an anti-slip mat item no. 628050-01, 50x30 cm is included, weight 2 kg.

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For all exercises please stand relaxed, upright, with lightly bended knees. Take care about a correct foot position: (fig. b) Fore foot turns inwards – weight is on big toe joint. Back foot (heel) turns outwards. Now adjust your leg-axle by turning the lower leg inwards and the thigh outwards to the hip resulting in a straight line of ankle joint – knee joint – hip joint (fig. a). Please always make sure to have a correct body posture.

Stability comes out of the feet. All foot and leg muscles along the leg axle need to collude. The more safer and well controlled your stand is, the more stable your trunk, the merrier your body posture and movement. While doing the exercises you will see there is a holistic effect.

Important! Make sure that exercises are always done controlled, slowly and even fluently

Variants during exercises:

- Change your head position: Turn your head left and right, look up and down.
- with opened or closed eyes.
- Balance a ball, beanbags or a table tennis racket additionally

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Exercises with the Footboards moveble fulcrum

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EXERCISES WITH THE FLAT SIDE UP:

Initial position:

Standing in lengthwise position on the Footboards.

Exercise scope:

- a) 10 sec. two-legged
- b) Seesaw back and forth for 10 sec. with both legs (fore foot/heel)
- c) Same as a) but each 10 sec. one-legged left/right

English

Initial position:

Position your feet such that the toes stand in front. Now tuck up your toes forcefully against the resistance of the Footboards.

Exercise scope:

- a) 5x in parallel
- b) 5x changing left/right
- c) 5x in parallel, in doing so lift up your heels

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Exercises with the Footboards moveble fulcrum

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English

Initial position:

Standing two-legged in cross position on one Footboard. Be aware of stretching your horizontally foot arch. Your foot will border the bulge of the Footboards more and more.

Exercise scope:

- Change the position of your feet several times for wards/backwards. Let each position sink in for 3 sec.
- Lift up you legs alternately. Tread on your feet.
- Build up tension in your feet and pull your toes and heel upwards for 5 times.

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Exercises with the Footboards moveble fulcrum

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English

Initial position:

Standing two-legged in cross position on the first Footboard. Now change slowly, smooth and controlled to the second Footboard and back. Let each position sink in for 3 sec.


Exercise scope:

- Left forwards, right follows, left backwards, right follows.
- Right forwards, left follows, right backwards, left follows.

a)

a)

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Exercises with the Footboards moveble fulcrum

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EXERCISES WITH THE FLAT SIDE UP:

Initial position:

Standing in lengthwise position on the Footboards.

Exercise scope:

- 10 sec two-legged, straightened, stable
- 5 x seesawing inwards (shift your weight on the inner side of your foot) and outwards (shift your weight on the outer side of your foot)
- 5 x seesawing two-legged in parallel left/right (skiing)
- Same as a) but each 10 sec one-legged left/right

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Exercises with the Footboards moveble fulcrum

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
English

Initial position:

Standing two-legged in cross position on one Footboard.

Exercise scope:

- 5x lift up your heels to a horizontal position, keep it for 3 sec and lower it again.
- 5x stand up from horizontal heel position into toe stand, keep it for 3 sec and lower in horizontal position. Don't set down your heels in between times.

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Exercises with the Footboards moveble fulcrum

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
English

Initial position:

Standing two-legged in cross position on one Footboard. Lift up your heels to a horizontal position. Now...

Exercise scope:

- a) ... 5x changing in an one-legged-stand left/right. Let each position sink in for 3 sec. Don't set down your heels.
- b) ...5x changing in a toe stand left/right. Let each position sink in for 3 sec. Don't set down your heels.

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Exercises with the Footboards moveble fulcrum

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
English

Initial position:

Standing two-legged in cross position on the first Footboard. Lift up your heels. Now change slowly, smooth and controlled to the second Footboard and back. Let each position sink in for 3 sec.

Exercise scope:

- a) Left forwards, right follows, left backwards, right follows.
- b) Right forwards, left follows, right backwards, left follows.

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Exercises with the Footboards moveble fulcrum

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English

Initial position:

Standing two-legged in cross position on one Footboard.

Exercise scope:

- 5x lift up your fore foot to horizontal position, keep it for 3 sec and lower again.
- 5x change from horizontal fore foot position into heel stand (lift up your fore foot), keep it for 3 sec and lower back in horizontal position. Don't set down your toes in between times.

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Exercises with the Footboards moveble fulcrum

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
English

Initial position:

Standing with crossed legs lengthwise on the Footboards.

Exercise scope:

- Now shift your weight on the front crossed leg. Let your back crossed leg follow around your supporting leg in front, to find position again on the Footboard. Repeat for 4 times.
- Same as a) but now with changed initial position. The previous front leg starts now from back position.

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Exercises with the Footboards moveble fulcrum

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English

Initial position:

Put both Footboards lengthwise one after another. Place your foot centric on both Footboards.

Exercise scope:

- In one-leg-stand: Try to turn your foot in the metatarsus area, thus to mobilise it by turning the fore foot inwards and the heel outwards. Mobilise each foot for 15 sec.
- Standing with both feet one after the other on the Footboards (rope dance) that way, that each foot is placed on one Footboard. Change the foot position. Balance each foot position for 15 sec. (right foot/left foot in front).

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Further information on the issue of „A healthy foot“ on our website www.pedalo.de or in the brochure with the topic „A healthy foot“.

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Dear customer,

The pedalo®-article you have chosen is made of natural wood. Differences in colour, structure, grain, small included knobs or scars prove the genuineness of our product. They are natural and therefore features of using original wood and no quality reduction.

Surface:

Delivered untreated. If you wish you can oil the solid ash wood Footboards moveable fulcrum with every natural oil or just leave it raw – as is. To have the Footboards moveable fulcrum slightly griped you might rub the surface with a little water one time, so the wooden fibers straighten up and the surface will have more grip.

Please note:

For the use of this device a certain degree in motor abilities and circumspection is assumed. In case of doubt first use device only with assistance.

copyright by Holz-Hoerz GmbH, printed copies also partially are not allowed.

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