SISSEL® Myofascia Roller

- ▮ Gebrauchsanweisung
- ▮ Instructions for use
- ▮ Mode d'emploi
- ▮ Gebruiksaanwijzing
- ▮ Instrucciones de uso
- ▮ Istruzioni per l'uso
- ▮ Bruksanvisning
- ▮ 사용방법

SISSEL® Myofascia Roller

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SISSEL® Myofascia Roller

Treat muscle tensions, promote regeneration –

Effective, safe, skin-friendly, mobile

Perfect for an effective self-massage as well as for professional use, loosens hardened muscles and sticky fascia, especially of the lower extremities (lower leg and thigh musculature). Regular use prevents continuing muscle tightening. The oval recess is ergonomic and reduces pressure on the wrists while exercising.

Extremely stable, for years of use

Generously sized (40 x 15 cm) for sufficient pelvis support

Very effective in combination with the SISSEL® Vitalyzer or SISSEL® Spiky Ball.

Recommended use

- Lower back pain.
- Movement impairments
- Diffuse pain in knee, hip or shoulder after consultation with the treating medical professional
- Application supports the correction of motion sequences and incorrect posture and can be used for balance training.

- Used as therapeutic aid in physio therapy, osteopathy, massage of connective tissue and sports medicine.

Contraindications

- Injuries, fractures, open wounds
- Osteoporosis
- Varicose veins

To perform effectively and free of pain, our musculature or actually the surrounding fascial tissue must be loosed specifically. Fascia (fibrous, connective tissue) knotted by injury or degeneration is noticeably loosened with the SISSEL® Myofascia Roller, a prerequisite for persistent muscle release. Afterwards the musculature is ready for stretching, loosening exercises and more workout. The active loosening of tensions throughout the body, helps target back pain in particular. Made in Germany from high grade EPP, recyclable, hygienic, easily cleaned with a moist towel.

Registered medical device in accordance with EU directive 93/42

Note:

Initially, and depending on the degree of tenseness the effect of the SISSEL® Myofascia Roller can be intense. We recommend using a cushioning exercise mat. Stop using if in pain and consult your doctor.

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SISEL® GmbH