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DEUTSCH

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
+49 8581 9699340
hello@triggerbow.com

ENGLISH

THANK YOU VERY MUCH FOR CHOOSING OUR TRIGGERBOW.

You are now holding an excellent quality product in your hands, with which you can immediately start a pain-free life. Before you start with the self-treatment, please read this manual carefully and also keep it for future reference. Please make this manual accessible to all other users of your triggerbow and pay attention to the following instructions.

Safety instructions

Please use the triggerbow for self-massage only. The product is intended for private use only and the treatment method specified in the instructions for use must be followed. The 22 mm stainless steel ball serves as a treatment piece. It is held in the adapter device by a magnet. This feature ensures a comfortable rolling function, but the ball may fall out of the holder in case of lateral or excessive force.

The manufacturer is not liable for any damage caused by such improper use.

Do not use the triggerbow ...

- immediately after surgery,
- during pregnancy,
- in case of injury or infection in the body region to be massaged (e.g. open wounds, spondylarthrosis),
- in case of pain of unexplained cause,
- in case of impaired perception (e.g. pain-relieving medications, alcohol),
- in case of severe diseases,
- in case of tumors in the treatment area,
- when taking anticoagulant medication,
- immediately after operations in the spinal area,
- in case of severe nerve diseases or acute inflammations.

WARNING

IF YOU ARE UNSURE WHETHER YOU CAN TREAT YOURSELF WITH THE TRIGGERBOW IT IS BEST TO CONSULT YOUR DOCTOR.

SAFETY INSTRUCTIONS

Please note:

- Patients with pacemakers or implantable cardioverter defibrillators must not position the treatment center in the area of the heart.
- The packaging material must be kept away from children.
- Repair the device only if proper functioning is guaranteed afterwards.
- The triggerbow can be used from the age of 16.
- Use from middle childhood (completion of 7th year of life) must be done under the supervision of an adult.
- Caution: Danger of suffocation due to small parts (stainless steel ball) or packaging material.
- Please use the triggerbow only if it is safe to use, i.e. with correct application and functioning device.

Notes on the integrated magnet

Normally, your body will not come into contact with the integrated magnet when using the triggerbow, as it is securely glued. However, we would like to give you some general warnings about the use of magnets to avoid possible injuries.

- **Electronic devices:** Hard disks and EC cards as well as other electronic devices such as TVs, monitors and watches are susceptible to strong magnetic fields. Loss of data and damage to the devices can be the result of contact with a magnet. Therefore it is important to keep a safe distance from these devices (at least 1 meter). In addition, magnets must not be placed near meters such as gas meters.

SAFETY INSTRUCTIONS

- **Pacemakers:** Magnets can also affect the function of pacemakers, so a safe distance is absolutely necessary.
- **Children:** Magnets do not belong in children's hands, for example, there is a risk of swallowing. So make sure to store the magnets in places that are not accessible to children.
- **Processing:** Magnets must not be sawn or drilled through. They must also not be machined with files. The magnets could break and the dust that would be created during processing is very flammable. Also keep away from open flames, as magnets can lose their power above a certain temperature.
- **Risk of crushing:** Magnets require careful handling, because there is a high risk of crushing, from which even gloves will not protect you. Never try out magnets and their effects on parts of people's bodies, e.g. on their ears.
- **Shattering hazard:** Magnets can also shatter when they hit the ground or collide with each other. Wear safety gloves and goggles for your protection.
- **Allergies:** Many magnets are coated with nickel. Constant contact with nickel can cause a nickel allergy. If you have such an allergy, you should avoid touching magnets or wear gloves.

It is essential: Magnets should never be left lying around carelessly. In the worst case, uninformed people can injure themselves because they are not aware of the magnetic force and its dangers. If you are giving magnets or a device with an integrated magnet as a gift or as a gift to someone else, please include these warnings and point out the possible dangers!


triggerbow disclaims any liability for damages caused by improper use of the magnets.

Application

Muscle stiffness occurs quickly and is usually a very painful issue. It's great that you can now treat your muscle stiffness yourself! The best way to do this is to free your upper body and apply your triggerbow directly to the skin. This is how you proceed:

- By independently **palpating** the affected muscles, you will find the areas that are slightly to severely painful and show increased tension.
- At these points (trigger points), for example, you perform the **pressure technique**. Use only enough force so that the perceived stimulus still has a pleasant character.
- Always maneuver the massager with **both hands** so that the treatment piece reaches the desired section of the back.
- You can also treat your back muscles with the triggerbow by a pleasant **rolling massage**. This technique is especially beneficial and relaxing for larger areas.

The device offers both a punctual and a flat massage with **individual pressure**. Increased muscle tension can thus be reduced. Stop the massage if you find it painful or uncomfortable.

Pressure technique

Press on the localized muscle stiffness and maintain this pressure for 40-60 seconds. You can repeat this maximum times.

Rolling movements

Massage the affected muscle area with circular rolling movements for a maximum of 10 minutes.

Application areas

GREEN

Treatment is here


RED

Do not apply pressure/
massage treatment to
firm, bony areas!

Further information on both techniques as well as specific application examples can be found online in our blog at triggerbow.com/blogs/news as well as on our therapy platform

 **triggerbow**^{PLUS}

plus.triggerbow.com



APPLICATION

Massage positions

The treatment is performed **sitting**, **standing** or **leaning against a wall**. Breathe calmly and try to relax. During distracting factors should not be present during the treatment.

You may massage yourself several times a day with the rolling technique, but please use the pressure technique only once a day for each trigger point. After 2-3 days you can already notice the first improvements. Reflect your results and adjust the treatment intensity accordingly.

APPLICATION

The picture shows the **wrong area of application!** With the triggerbow you treat only the soft tissues (muscles) of the human body. Pressure treatment should not be applied to **sensitive and bony parts of the body**, e.g. the spine (spinous processes), iliac crest, sacrum, tibia, cranial bones or ribs.

In addition, please pay attention to an **adapted intensity**, otherwise unwanted side effects may occur, e.g. bruises (hematomas).

The triggerbow does not replace therapeutic treatment. If you have any questions, please contact us or your doctor.

do not use in
the bone area

Online therapy platform

We are happy to support you in the self-treatment of your trigger points with detailed step-by-step instructions for each muscle area. No matter what kind of pain you suffer from - from headache to neck pain to back pain: On our **therapy platform triggerbow® PLUS** you will find the right video tutorial and much more information. Just have a look!

plus.triggerbow.com


Guides

So much for the basic application of the triggerbow. We would like to show you all the pain areas of your body in detail, including the trigger points that occur there as well. With the help of our detailed **guides** and the body maps they contain, you can then treat the following pains even more specifically:

- **headaches,**
- **neck and shoulder pain,**
- **back pain and**
- **pain in the foot**

A lot of additional information and tips will make you a trigger point professional and master of self-treatment.

Just click through to our store at triggerbow.com and try it out!

Cleaning and maintenance

- Hot and wet can damage the wood surface.
- Do not use chemical cleaners.
- Clean the treatment adapter with a damp cloth and then dry it.
- Do not store the unit in a damp environment, but preferably at room temperature.
- Any contact with damp surfaces may cause unwanted changes in the shape of the product.
- The product is made of wood and therefore easily combustible.
- The life of the triggerbow depends on careful handling.
- Do not expose the device to strong temperature fluctuations, shocks or chemicals.
- Do not use abrasive cleaning agents.
- The cleaning process should be carried out by an adult.

Disposal

The device should not be disposed of in household waste at the end of its service life due to the integrated magnet. Disposal can be carried out via appropriate collection points.

Contact

If you have any questions or problems, please feel free to contact us:

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vivadesign GmbH