## aero step® XL & aero step® XL sensitive

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### Getting started - finding your balance

Standing on the Aero-Step® XL, intensify the upright posture and tense the body lengthwise, with intentional awareness of the cushion. Then move the feet forwards, backwards and to the sides, followed by circling the knees or pelvis in both directions. Then remain standing quietly, not tipping over forwards or backwards, finding your balance, perhaps closing your eyes and becoming aware of the act of balancing.

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### Criss-cross exercise - linking and more

Standing on the Aero-Step® XL, firmly touch the opposing leg with alternating hands. For greater intensity, accompany the exercise with a movement of the head to the left and right.

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### Overall body tensing

Standing on the knees with an upright, neutral posture. Tip upper body forwards and backwards, holding it in a stable posture. To raise the level of difficulty, lift the feet.

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### Back - buttocks

Lie down in a comfortable position, face down, on the Aero-Step® XL. Stretch out legs until they are raised from the floor, then raise the stretched upper body and maintain this lengthwise tension. In this position, raise and lower the legs alternately and/or rotate the upper body slightly to the left and right.

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### Shoulders - arms

On all fours, support the lower arms on the Aero-Step® XL. Shifting the weight to one side, raise the opposing arm laterally to the side and diagonally forwards.

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### Relaxation - stretching

A pleasant and relaxing stretching session is recommended after the strong stimuli of the exercises: Allow the breath to ebb and flow.

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GB

The exercises are selected according to functional criteria for your health and wellbeing.

#### Training tips:

- Always begin by assuming the initial position very precisely and tensing the lower abdomen somewhat
- Perform all of the exercises slowly
- Increase the number of repetitions to 20 per exercise
- Breathe quietly and naturally during the training session
- Leave out exercises that are unpleasant
- Rest between exercises
- »Wagging on the cushion is something you should look forward to: it will hasten your progress

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